

Product datasheet for **SC123964**

KCNA5 (NM_002234) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KCNA5 (NM_002234) Human Untagged Clone
Tag:	Tag Free
Symbol:	KCNA5
Synonyms:	ATFB7; HCK1; HK2; HPCN1; KV1.5; PCN1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_002234, the custom clone sequence may differ by one or more nucleotides

```

ATGGAGATCGCCCTGGTGCCCTGGAGAACGGCGGTGCCATGACCGTCAGAGGAGGCGATGAGGCCCGGG
CAGGCTGCGGCCAGGCCACAGGGGAGAGCTCCAGTGTCCTCCCGACGGCTGGGCTCAGCGATGGGCCCAA
GGAGCCGGCGCCAAAGGGGCGCGCGCAGAGAGACCGGACTCGGGAGTGCGGCCCTTGCTCCGCTG
CCGACCCCGGGAGTGCGGCCCTTGCTCCGCTGCCAGAGGAGCTGCCACGGCCTCGACGGCCGCTCCCG
AGGACGAGGAGGAAGAAGGCGATCCCGGCTGGGCACGGTGGAGGACCAGGCTCTGGGCACGGCGTCCCT
GCACCACCAGCGCTCCACATCAACATCTCCGGGCTGCGCTTTGAGACGCAGCTGGGCACCCTGGCGCAG
TTCCCAAACACTCTGGGGGACCCCGCAAGCGCTGCGCTACTTCGACCCCTGAGGAACGAGTACT
TCTTCGACCGCAACCGGCCAGCTTCGACGGTATCCTCTACTACTACCAGTCCGGGGGCCGCTGCGGAG
GCCGTC AACGTCTCCCTGGACGTGTCGCGGACGAGATACGCTTCTACCAGTGGGGGACGAGGCCATG
GAGCGCTTCGCGAGGATGAGGGCTTCATTAAGAAGAGGAGAAGCCCTGCCCGCAACGAGTTCAGC
GCCAGGTGGCTTATCTTCGAGTATCCGAGAGCTCTGGTCCGCGCGGCCATCGCCATCGTCTCGGT
CTTGTTATCCTCATCTCCATCATCACCTTCTGCTTGGAGACCCTGCCTGAGTTCAGGGATGAACGTGAG
CTGCTCCGCCACCTCCGGCGCCACCAGCCTCCCGCGCCCGCCCTGGGGCAACGGCAGCGGGTCA
TGGCCCGCCCTCTGGCCCTACGGTGGCACCGCTCCTGCCAGGACCCTGGCCGACCCCTTCTTCATCGT
GGAGACCAGTGCCTCATCTGGTTCACCTTCGAGCTGCTCGTGCCTTCTTCGCTGCCCGCAAGGCA
GGGTTCTCCCGAACATCATGAACATCATCGATGTGGTGGCCATCTCCCTACTTCATCACCTGGGCA
CCGAACGGCAGAGCAGCAGCCAGGGGTGGAGGAGGCGCCAGAATGGGCAGCAGGCCATGTCCCTGGC
CATCTCCGAGTATCCGCTGGTCCGGGTGTTCCGCATCTTCAAGCTCTCCCGCACTCCAAGGGGCTG
CAGATCTGGCAAGACCTTCAGGCTCCATGAGGGAGCTGGGGCTGCTCATCTTCTTCTTCTTCTG
GGTTCATCTTCTCCAGTGCCGTCTACTTCGCAGAGGCTGACAACCAGGGAACCCATTCTCTAGCAT
CCCTGACGCTTCTGGTGGGAGTGGTACCATGACCACTGTGGGCTACGGGGACATGAGGCCCATCACT
GTTGGGGCAAGATCGTGGGCTCGTGTGTGCCATCGCCGGGTCCTCACCATTGCCCTGCCTGTGCCCG
TCATCGTCTCAACTTCAACTACTTCTACCACCGGAAACGGATCACGAGGAGCCGGCAGTCTTAAAGGA
AGAGCAGGGCACTCAGAGCCAGGGGCGGGGCTGGACAGAGGAGTCCAGCGGAAGTCCAGCGGAGCAGG
GGATCTTCTGCAAGGCTGGGGGACCTGGAGAATGCAGACAGTCCCGAAGGGGCAGTGCCTCCCTAG
AGAAGTGAACGTCAAGGCAAGAGCAACGTGGACTTGGGAGGTCCTTTATGCCCTCTGCCTGGACAC
CAGCCGGAAACAGATTTGTGA
    
```

5' Read Nucleotide Sequence:

>OriGene 5' read for NM_002234 unedited

```

NTTTCCTGGAGAACGGGTGCCATGACCGTCAGAGGAGGCGATGAGGCCCGGGCAGGCT
GCGGCCAGGCCACAGGGGAGAGCTCCAGTGTCCTCCCGACGGCTGGGCTCAGCGATGGGC
CCAAGGAGCCGGCGCAAAGGGGCGCGCGCAGAGAGACCGGACTCGGGAGTGCGGC
CCTTGCTCCGCTGCCGACCCGGGAGTGCGGCCCTTGCTCCGCTGCCAGAGGAGTGC
CACGGCTCGACGGCCGCTCCCGAGGACGAGGAGGAAGAAGGCGATCCCGGCTGGGCA
CGGTGGAGGACCAGGCTCTGGGCACGGCTCCCTGCACCACCAGCGCTCCACATCAACA
TCTCCGGGCTGCGCTTTGAGACGCAGCTGGGCACCCTGGCGAGTTCCTCCAAACACTCC
TGGGGGACCCCGCAAGCGCTGCGCTACTTCGACCCCTGAGGAACGAGTACTTCTTCG
ACCGCAACCGGCCAGCTTCGACGGTATCCTCTACTACTACCAGTCCGGGGGCCGCTGC
GGAGCCGGTCAACGTCTCCCTGGACGTGTTCCGGGACGAGATACGCTTCTACCAGCTGG
GGGACGAGGCCATGGAGCGCTTCGCGAGGATGAGGGCTTCATTAAGAAGAGGAGAAGC
CCCTGCCCGCAACGAGTTCAGCGCCAGGTGTGGCTTATCNTCGAGTATCCGGAGAGCT
CTGGGTCCGCGCGGCCATCGCCATCGTCTCGGTCTTGTTATCTCATCTCCTCATCACC
TTCTGCTTGAGACCNTGCTGAGTTAGGGATGACGTGAGCTGCTCCN
    
```

3' Read Nucleotide Sequence:	>OriGene 3' read for NM_002234 unedited NTTGCCTGGATGGCACTTCCAGGCCAGGAGAGGCACTGGGGAGGGGTCACAGGGATGCCA CCCGGGATCTGTTTCAGGAAACAGCTATGACCGCGGCCCAATCTAGAGTCGAGTTTTCTT TCTTTGGGTTTTTTTGAACAAGGTAACGGTTTTATTTCTGGGATAAAAAATAAGACAAT AAATTATCAGAAGTATTCCTCCAATATCATTAAATGCTTTTTGAAGAAAGTTACAAGAA TTGTTTAAAGGCGCACAGAAAAAGTCAGATTTGCCTTTCTTCAAGAACATACTCTGCCCTA GACACAAAAGCACCAATACCACATCCCAGAAAAGACAGAGGCTTGGAGACACATGAAACAG ATCAGCCCATACCAAGATCTGACTCTAATATAAGTAACACAGAAAAACATCGGGACAA CACTTTCTCCATAAAATACAACATGTATCTTGAAAATATACCCAAAATATCTCTATCTGA ACTACAGTACAAAAAGCCACTTACCAGGCATAAAAAATAATTTTTCCCCTATACACGAAT ACAATCATATGTCAAGCAAATTTGGAACACATTGATGGGAGAGTTACATAATCTCTTTGG ATCTGGGCTATTTAGAGTCATTTTTGTTGGGGTATGCCTCTGCTTGGGTATACGAGACT GGCTACAAATACATNCATTCTAAAGCGGTTAACAATTTCAACAACCTCTCGGTTATTGG AATGGGGTACTTTTTAAAAATAAGGCACACCTTGGTGTGTGGAAAGGAATTATCAAATT GGGGGGCTCTCGGGGGGGGGAAAAAGGTTTAAAGGGAGGGGTATCTGNG
Restriction Sites:	NotI-NotI
ACCN:	NM_002234
Insert Size:	3000 bp
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.
	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_002234.2 , NP_002225.2
RefSeq Size:	2865 bp
RefSeq ORF:	1842 bp

Locus ID:	3741
UniProt ID:	P22460
Cytogenetics:	12p13.32
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
Gene Summary:	<p>Potassium channels represent the most complex class of voltage-gated ion channels from both functional and structural standpoints. Their diverse functions include regulating neurotransmitter release, heart rate, insulin secretion, neuronal excitability, epithelial electrolyte transport, smooth muscle contraction, and cell volume. Four sequence-related potassium channel genes - shaker, shaw, shab, and shal - have been identified in <i>Drosophila</i>, and each has been shown to have human homolog(s). This gene encodes a member of the potassium channel, voltage-gated, shaker-related subfamily. This member contains six membrane-spanning domains with a shaker-type repeat in the fourth segment. It belongs to the delayed rectifier class, the function of which could restore the resting membrane potential of beta cells after depolarization and thereby contribute to the regulation of insulin secretion. This gene is intronless, and the gene is clustered with genes KCNA1 and KCNA6 on chromosome 12. Defects in this gene are a cause of familial atrial fibrillation type 7 (ATFB7). [provided by RefSeq, May 2012]</p>