

Product datasheet for **SC124621**

TAF6 (NM_139315) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAF6 (NM_139315) Human Untagged Clone
Tag:	Tag Free
Symbol:	TAF6
Synonyms:	ALYUS; MGC:8964; TAF(II)70; TAF(II)80; TAF2E; TAFII-70; TAFII-80; TAFII70; TAFII80; TAFII85
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_139315, the custom clone sequence may differ by one or more nucleotides

```
ATGGCTGAGGAGAAGAAGCTGAAGCTTAGCAACACTGTGCTGCCCTCGGAGTCCATGAAGTGGTGGCTG
AATCCATGGGCATCGCCAGATTAGGAGGAGACCTGCCAGCTGCTAACGGATGAGGTGAGTACCCGCAT
CAAAGAGATCGCACAGGATGCCTTGAAGTTCATGCACATGGGGAAGCGGAGAAGCTCACCACCAGTGAC
ATTGACTACGCCTTGAAGCTAAAGAATGTCGAGCCACTCTATGGCTTCCACGCCAGGAGTTCATTCCTT
TCCGCTTCGCCTCTGGTGGGGCCGGGAGCTTTACTTCTATGAGGAGAAGGAGGTTGATCTGAGCGACAT
CATCAATACCCCTCTGCCCGGGTGGCCCTGGACGTCTGCCTCAAAGCTCATTGGCTGAGCATCGAGGGC
TGCCAGCCAGCTATCCCCGAGAACCCGCCCCAGCTCCCAAAGAGCAACAGAAGGCTGAAGCCACAGAAC
CCCTGAAGTCAGCCAAGCCAGGCCAGGAGGAAGACGGACCCCTGAAGGGCAAAGGTCAAGGGGCCACCAC
AGCCGACGGCAAAGGGAAAGAGAAGAAGGCGCCGCCCTTGGTGGAGGGGGCCCCCTTGGACTGAAGCCC
CGGAGCATCCACGAGTTGTCTGTGGAGCAGCAGCTCTACTACAAGGAGATCACCGAGGCCCTGCGTGGCT
CCTGCGAGGCCAAGAGGGCGGAAGCCCTGCAAAGCATTGCCACGGACCCCTGGACTGTATCAGATGCTGCC
ACGGTTACGTACCTTTATCTCGGAGGGGTCCGTGTGAACGTGGTTCAGAAACACCTGGCCCTACTCATC
TACCTGATGCGTATGGTAAAGCGCTGATGGACAACCCACGCTCTATCTAGAAAAATACGTCCATGAGC
TGATTCAGCTGTGATGACCTGCATCGTGAGCAGACAGTTGTGCCTGCGACCAGATGTGGACAATCACTG
GGCACTCCGAGACTTTGCTGCCCGCTGGTGGCCAGATCTGCAAGCATTTTAGCACAACCACTAACAAC
ATCCAGTCCCGGATCACAAGACCTTACCAAGAGCTGGGTGGACGAGAAGACGCCCTGGACGACTCGTT
ATGGCTCCATCGCAGGCTTGGCTGAGCTGGGACACGATGTTATCAAGACTCTGATTCTGCCCCGGCTGCA
GCAGGAAGGGGAGCGGATCCGCAGTGTGCTGGACGGCCCTGTGCTGAGCAACATTGACCGGATTGGAGCA
GACCATGTGCAGAGCCTCTGCTGAAACACTGTGCTCCTGTTCTGGCAAAGCTGCGCCACCCGCTGACA
ATCAGGACGCCTATCGGGCAGAATTCGGGTCCCTTGGGCCCTCCTCTGCTCCCAGGTGGTCAAGGCTCG
GGCCAGGCTGCTCTGCAGGCTCAGCAGGTCAACAGGACCACTCTGACCATCACGCAGCCCCGGCCACG
CTGACCCCTCTCGCAGGCCCCACAGCCTGGCCCTCGCACCCCTGGCTTGTGAAGTTCTGGCTCCATCG
CACTTCTGTCCAGACTGGTGTCTGCACGAGCGGCTGCCCCACCACAGCCTTCCCCTCTCAACCAA
GTTTATTGTAATGTCATCGTCTCCAGCGCCCATCCACCAGCAGGTCTGTCCCTCAGCACCTCGGCC
CCCGGCTCAGGTTCCACCACCACTTCGCCCGTACCACCACCGTCCCAGCGTGCAGCCCATCGTCAAGT
TGGTCTCCACCGCCACCACCGCACCCCCAGCACTGCTCCCTCTGGTCTGGGAGTGTCCAGAAGTACAT
CGTGGTCTCACTTCCCCAACAGGGGAGGGCAAAGGAGGCCCCACCTCCCATCCTTCTCCAGTTCCTCCC
CCGGCATCGTCCCGTCCCCTCAGCGGCAAGTGCCTTTGTGGGGGAAGCAGGAGGCTGGGGACAGTC
CCCTCCAGCTCCAGGACTCCAAAAGCAATGGCTCCCAGCCAACTCCGGCTCCCCTCAGCCTGCTCC
GTGA
```

5' Read Nucleotide Sequence:	>OriGene 5' read for NM_139315 unedited GAATACGACTCACTATAGGGCGGCCGGAATTCGGCACGAGGGGACTCCAGAATGGCTGA GGAGAAGAAGCTGAAGCTTAGCAACACTGTGCTGCCCTCGGAGTCCATGAAGGTGGTGGC TGAATCCATGGGCATCGCCAGATTTCAGGAGGAGACCTGCCAGCTGCTAACGGATGAGGT CAGCTACCGCATCAAAGAGATCGCACAGGATGCCTTGAAGTTCATGCACATGGGGAAGCG GCAGAAGCTCACCACCAGTGACATTGACTACGCCTTGAAGCTAAAGAATGTCGAGCCACT CTATGGCTTCCACGCCAGGAGTTCATTCTTTCCGCTTCGCCTCTGGTGGGGGCCGGA GCTTTACTTCTATGAGGAGAAGGAGTTGATCTGAGCGACATCATCAATACCCCTCTGCC CCGGGTGCCCTGGACGTCTGCCTCAAAGCTCATTGGCTGAGCATCGAGGGCTGCCAGCC AGCTATCCCCGAGAACCCGCCCCAGCTCCCAAAGAGCAACAGAAGGCTGAAGCCACAGA ACCCCTGAAGTCAGCCAAGCCAGGCCAGGAGGAAGACGGACCCCTGAAGGGCAAAGGTCA AGGGGCCACCACAGCCGACGGCAAAGGGAAAGAGAAGAAGGCCGCCCTTGCTGGAGGN NGGCCCCCTTGCAGTGAAGCCCGAGCTCCACGAGTNGTCTGTGGAGCAGCAGCTCTAC TACAAGGAGATCACCGAGCCTGCGTGGGCTCCTGCGAGCCAAGAGGGCGGAGCCCTGCAA GCATTGCACGGACCCTGGACTGTATCAGAGCTGNCCCGNTTCAGTACTTATCTCGGAGGG GGTGCTGGGAACGTGTTCAAACAACCTGGCTACTCATTACTTGATGCTTGGGNGAAGCCT GATGACACCCCGCTCTACAGNAAAACTTCTC
Restriction Sites:	NotI-NotI
ACCN:	NM_139315
Insert Size:	2500 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_139315.1 , NP_647476.1
RefSeq Size:	2489 bp
RefSeq ORF:	2034 bp
Locus ID:	6878
UniProt ID:	P49848
Cytogenetics:	7q22.1
Protein Families:	Transcription Factors

Protein Pathways: Basal transcription factors

Gene Summary: Initiation of transcription by RNA polymerase II requires the activities of more than 70 polypeptides. The protein that coordinates these activities is transcription factor IID (TFIID), which binds to the core promoter to position the polymerase properly, serves as the scaffold for assembly of the remainder of the transcription complex, and acts as a channel for regulatory signals. TFIID is composed of the TATA-binding protein (TBP) and a group of evolutionarily conserved proteins known as TBP-associated factors or TAFs. TAFs may participate in basal transcription, serve as coactivators, function in promoter recognition or modify general transcription factors (GTFs) to facilitate complex assembly and transcription initiation. This gene encodes one of the smaller subunits of TFIID that binds weakly to TBP but strongly to TAF1, the largest subunit of TFIID. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2010]

Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Both variants 1 and 2 encode the same isoform (alpha).