

Product datasheet for **SC124743**

GOLGA6L9 (NM_198181) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA6L9 (NM_198181) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA6L9
Synonyms:	GOLGA6L20
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_198181, the custom clone sequence may differ by one or more nucleotides

ATGTGGCCCCAACCCCGCTCCCTCCCCACCCCGCATGTTCAGAAAAACACAACAGGGGAAATTGGCCG
CAGCCAAGAAAAAGTTAAAGCATATTGGCAGAGGAAGAGCCCTGGCATTCCAGCAGGAGCTAACAGGAA
AAAGAAAGTCAATGGCAGTAGCCCTGACACAGCCACTTCTGGTGGTTACCACTCACCTGGGGATTACGCA
ACAGGTATCTACGGGGAGGGCCGTGCATCCTCTACTACCCTGCAGGATCTGGAGAGCCAGTACCAAGAAC
TAGCAGTAGCCCTGGATTCAAGCTCCGCAATAATCAGTCAACTCACTGAAAAATCAATTCACTGGTTTCG
CACATCTAAGGAGGAGAAGAAGCATGAGATACATCTGGTACAGAAGCTTGGGAGGAGCTTGTTCAAACTC
AAAAACCAGACGGCTGAACCCCTGGCCCCACAGCCCCCAGCAGGGCCATCTAAGGTGGAGCAGCTACAAG
ATGAGACCAACCCTAAGGAAGGAGCTAGAGAGTGTGGGAAGACAGCTCCAGGCTGAGGTGGAAAAACA
TCAGATGTTGAGTCTCCTGAACAGGAGACAGGAGGAGAGGCTACGTGAACAGGAGGAGAGGCTACGTGAA
CAGGAGGAGAGGCTGTGTGAACAGGAGGAGAGGCTGTGTGAACAGGAGGAGAGGCTACGTGAACAGGAGG
AGAGGCTGTGTGAACAGGAGAAGCTGCCAGGGCAGGAGAGGCTGCTGGAAGAGGTGGAGAAGCTGTTAAA
ACAGGAGAGGGCGCAGGAGGAGCAGGAGAGGCTGCTGGAGAGGGAGAGGCTGCTGGACGAGGTGGAGGAG
CTCCTGGAGCAGGAGAGGCTTCGGCAACAGGATGAGAGGCTGTGGCAGCAGGAGACTCTGCGGGAGCTGG
AGAGGTTGCGGGAGCTGGAGAGGATGCTGGAGCTGGGGTGGGAAGCCCTGTACGAGCAGCGGGCCGAGCC
ACGCAGCGGCTTCGAGGAGCTGAACAACGAGAACAAGAGCACACTGCAGTTGGAGCAGCAAGTAAAGGAG
CTGGAGAAGTCGGGTGGAGCTGAAGAGCCAAGAGGCTCCGAGTCTGCAGCAGCAGCCAGACCAGTACCTG
GAGCCCCAGTCCCAAGGAGCTTGATGTGCGGACAAGCAGGCTGGACTCCCCAGGAGCAGCCAGGCTT
GAGTGGAGAAGCTGTTGGTACAGGAGAGGCGGCAGGAGGAGCAGGAGAGGCTGCATGCCATTCTTTTCG
GCTGCCGAGAACAGGAGCTAACATCACCATCATCTAA



[View online »](#)

5' Read Nucleotide Sequence:

[illegible]

Restriction Sites:

Notl-Notl

ACCN:

NM_198181

OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:

NM 198181.1, NP 937824.1

RefSeq Size:

4375 bp

RefSeq ORF:

825 bp

Locus ID:

440295

UniProt ID:

A6NEM1

Cytogenetics:

15q25.2