

Product datasheet for **SC124776**

MRPL21 (NM_181515) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MRPL21 (NM_181515) Human Untagged Clone
Tag:	Tag Free
Symbol:	MRPL21
Synonyms:	L21mt; MRP-L21
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_181515, the custom clone sequence may differ by one or more nucleotides

ATGATCGTCACGGGGCAGTATGGCAGGCTCTTTGCCGTGGTGCACCTTTGCCAGCCGCCAGTGGAAGGTGA
 CCTCTGAAGACCTGATCTTAATTGGAAATGAACTAGACCTTGCGTGTGGAGAGAGAATTCGACTGGAGAA
 GGTCTGCTGGTTGGGGCAGACAACTTCACGCTGCTTGGCAAGCCACTCCTCGGAAAGGATCTTGTTCGA
 GTAGAAGCCACAGTCATTGAAAAGACAGAATCATGGCCAAGAATCATTATGAGATTGAGAAAAGGAAAA
 ACTTCAAGAAGAAAAGAATCGTCACGACCCCGCAGACTGTCCTCCGATAAACAGCATTGAGATTGCTCC
 GTGTTTGTGTGA

5' Read Nucleotide Sequence:

>OriGene 5' read for NM_181515 unedited
 GCACGAGCCTCGTGCCGAATTCGGCAGGAGGCGAGCATCTTCCCTGACGGTCACCTTAGG
 GCGGCTGGCGTCCGCGTGCAGCCACAGCATCCTGAGACCTTCGGGGCCCGGAGCAGCCTC
 CCTTTGGTCTGCTTCTCGAAGGTTCAATTCACAGAGCACTTCATATCTACCAGGATATGT
 TCCTAAACATCCCTGAGTTCACCACCTTGCCAGAAAGTTGTTCTGCCAGACCCAGTTGA
 GGAGACCAGACACCATGCAGAGGTCGTGAAGAAGGTGAATGAGATGATCGTCACGGGGCA
 GTATGGCAGGCTCTTTGCCGTGGTGCACCTTTGCCAGCCGCCAGTGGAAGGTGACCTCTGA
 AGACCTGATCTTAATTGGAAATGAACTAGACCTTGCGTGTGGAGAGAGAATTCGACTGGA
 GAAGGTCTGCTGGTTGGGGCAGACAACTTCACGCTGCTTGGCAAGCCACTCCTCGGAAA
 GGATCTTGTTGAGTAGAAGCCACAGTCATTGAAAAGACAGAATCATGGCCAAGAATCAT
 TATGAGATTGAGAAAAGGAAAACTTCAAGAAAAAAGAATCGTCACGACCCCGCAGACT
 GTCCTCCGATAAACAGATTGAA

Restriction Sites:	NotI-NotI
ACCN:	NM_181515


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_181515.1</u> , <u>NP_852616.1</u>
RefSeq Size:	726 bp
RefSeq ORF:	363 bp
Locus ID:	219927
UniProt ID:	<u>Q7Z2W9</u>
Cytogenetics:	11q13.3
Gene Summary:	<p>Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Multiple transcript variants encoding different isoforms were identified through sequence analysis although some may be subject to nonsense-mediated decay (NMD). [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) uses a different splice site in the 5' end, compared to variant 4, that likely leads to the use of a downstream start codon. The predicted protein (isoform a) is shorter than isoform d. The predicted ORF of this transcript has not been experimentally confirmed.</p>