

Product datasheet for SC125030

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

P2X5 (P2RX5) (NM_175081) Human Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: P2X5 (P2RX5) (NM_175081) Human Untagged Clone

Tag: Tag Free Symbol: P2X5

Synonyms: LRH-1; MGC47755; P2X5; P2X5R; purinergic receptor P2X, ligand-gated ion channel, 5;

purinergic receptor P2X5

Mammalian Cell

Selection:

None

Vector: pCMV6-XL4

E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >NCBI ORF sequence for NM_175081, the custom clone sequence may differ by one or more

nucleotides

ATGGGGCAGGCGGCTGCAAGGGGCTCTGCCTGTCGCTGTTCGACTACAAGACCGAGAAGTATGTCATCG CCAAGAACAAGAAGAAGTGGGCCTGCTGCAGCTGCTGCAGGCCTCCATCCTGGCGTACCTGGTCGTATG GGTGTTCCTGATAAAGAAGGGTTACCAAGACGTCGACACCTCCCTGCAGAGTGCTGTCATCACCAAAGTC AAGGGCGTGGCCTTCACCAACACCTCGGATCTTGGGCAGCGGATCTGGGATGTCGCCGACTACGTCATTC CAGCCCAGGGAGAACGTCTTTTTTTGTGGTCACCAACCTGATTGTGACCCCCAACCAGCGGCAGAACGT CTTGTGCTGAGAATGAAGGCATTCCTGATGGCGCGTGCTCCAAGGACAGCGACTGCCCACGCTGGGGAAGCG GTTACAGCTGGAAACGGAGTGAAGACCGGCCGCTGCTGCCGGAGAGAAACTTGGCCAGGGGCACCTGTG

AGATCTTTGCCTGGTGCCCGTTGGAGACAAGCTCCAGGCCGGAGTGA



5' Read Nucleotide Sequence: >OriGene 5' read for NM_175081 unedited

GTGTCAGGATTTTGTAATACGACTTCACTATAGGGCGGCCGCGAATTCGCACCAGGTCAA TAAAAGGCCGCCCGACCCGCCCCCGCGCCCCGCAGCCCTGCCGGACACCCCGGGCTGCA GCTGAGCGGGCGCAGACGGGCCGAGGCGCGCGCGCGCAGGGACCGAG TGCTCCCCATGAGCGCACGTGGGCCGGCGGCGCCAAGCCCGGCTGAGAGCGCGCCATG GGGCAGGCGGCTGCAAGGGGCTCTGCCTGTCGCTGTTCGACTACAAGACCGAGAAGTAT GTCATCGCCAAGAACAAGAAGGTGGGCCTGCTGTACCGGCTGCTGCAGGCCTCCATCCTG GCGTACCTGGTCGTATGGGTGTTCCTGATAAAGAAGGGTTACCAAGACGTCGACACCTCC CTGCAGAGTGCTGTCATCACCAAAGTCAAGGGCGTGGCCTTCACCAACACCTCGGATCTT GGGCAGCGGATCTGGGATGTCGCCGACTACGTCATTCCAGCCCAGGGAGAGAACGTCTTT TTTGTGGTCACCAACCTGATTGTGACCCCAACCAGCGGCAGAACGTCTGTGCTGAGAATG AAGGCATTCCTGATGGCGCGTGCTCCAAGGACAGCGACTGCCACGCTGGGGAAGCGGTTA CAGCTGGAAACGGAGTGAAGACCGGCCGCTGCCTGCGGAGAGAACTTGCCAGGGGCAC CTGTGAGATCTTTGCCTGGTGCCCGTTGGAGACAAGCTCCAGGCCGAAGGAGCATTTCCT GAAGGAGCCGGAGACTTCCCATTTTATTAAGACCCACATCGTTTCCCAAATTAACTTTTC AAAAGCATGTGATGGACGTCAAGAACAGATTTTTCTGAAATATGCCCTTTGGCCCAAGAC CCTACTGCCCATTTTCGATTGGTTC

Restriction Sites: Notl-Notl
ACCN: NM_175081

Insert Size: 2200 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a

point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative

RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 175081.1, NP 778256.1</u>

 RefSeq Size:
 958 bp

 RefSeq ORF:
 537 bp

 Locus ID:
 5026

 Cytogenetics:
 17p13.2

Protein Families: Druggable Genome, Ion Channels: ATP Receptors, Transmembrane

Protein Pathways: Calcium signaling pathway, Neuroactive ligand-receptor interaction





Gene Summary:

The product of this gene belongs to the family of purinoceptors for ATP. This receptor functions as a ligand-gated ion channel. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring downstream gene, TAX1BP3 (Tax1 binding protein 3). [provided by RefSeq, Mar 2011] Transcript Variant: This variant (3) lacks a segment in the coding region and differs in the 3' UTR, compared to variant 1. Its predicted isoform (C) has a shorter C-terminus, compared to isoform A.