

## Product datasheet for **SC125032**

### Uridine Phosphorylase 1 (UPP1) (NM\_181597) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Uridine Phosphorylase 1 (UPP1) (NM_181597) Human Untagged Clone
Tag:	Tag Free
Symbol:	Uridine Phosphorylase 1
Synonyms:	UDRPASE; UP; UPASE; UPP
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_181597, the custom clone sequence may differ by one or more nucleotides

```
ATGGCGGCCACGGGAGCCAATGCAGAGAAAGCTGAAAGTCACAATGATTGCCCCGTCAGACTTTTAAATC  
CAAACATAGCAAAAATGAAAGAAGATATTCTCTATCATTTCACCTACCAGTAGCAGACACAATTTCCC  
AGCCTTGTGGAGATGTGAAGTTGTGTGTGTTGGTGGAAAGCCCTCCCGGATGAAAGCCTTCATCAGG  
TGCGTTGGTGCAGAGCTGGGCCTTGACTGCCAGGTAGAGACTATCCCAACATCTGTGCGGGAAGTACC  
GCTATGCCATGTATAAAGTAGGACCGGTGCTGTCTGTCTGTCAGTCATGGTATGGGCATTCCTTCTATCTCAAT  
CATGTTGCATGAGCTCATAAAGCTGCTGTACTATGCCCGGTGCTCCAACGTCACATCATCCGATTGGC  
ACTTCTGGTGGGATAGGTCTGGAGCCCGGCACTGTGGTCATAACAGAGCAGGCAGTGGATACCTGCTTCA  
AGGCAGAGTTTGAGCAGATTGCTCTGGGGAAGCGGGTATCCGGAAAACGGACCTTAACAAGAAGCTGGT  
GCAGGAGCTGTTGCTGTGTTCTGCAGAGCTGAGCGAGTTCACCACAGTGGTGGGGAACACCATGTGCACC  
TTGGACTTCTATGAAGGGCAAGGCCGTCTGGATGGGGCTCTGCTCCTACACGGAGAAGGACAAGCAGG  
CGTATCTGGAGGCAGCCTATGCAGCCGGCGTCCGCAATATCGAGATGGAGTCTCGGTGTTTGCCGCCAT  
GTGCAGCGCCTGCGGCCTCCAAGCGGCCGTGGTGTGTGTACCCCTCCTGAACCGCCTGGAAGGGGACCAG  
ATCAGCAGCCCTCGCAATGTGCTCAGCGAGTACCAGCAGAGGCCCGCAGCGGCTGGTGTGAGCTACTTCATCA  
AGAAGAACTGAGCAAGGCCTGA
```

Restriction Sites:	NotI-NotI
ACCN:	NM_181597



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_181597.1</a> , <a href="#">NP_853628.1</a>
<b>RefSeq Size:</b>	1857 bp
<b>RefSeq ORF:</b>	933 bp
<b>Locus ID:</b>	7378
<b>Cytogenetics:</b>	7p12.3
<b>Protein Pathways:</b>	Drug metabolism - other enzymes, Metabolic pathways, Pyrimidine metabolism
<b>Gene Summary:</b>	<p>This gene encodes a uridine phosphorylase, an enzyme that catalyzes the reversible phosphorylation of uridine (or 2'- deoxyuridine) to uracil and ribose-1-phosphate (or deoxyribose-1-phosphate). The encoded enzyme functions in the degradation and salvage of pyrimidine ribonucleosides. [provided by RefSeq, Oct 2016]</p> <p>Transcript Variant: This variant (2) is the longer transcript. Both variants 1 and 2 encode the same isoform.</p>