

## Product datasheet for **SC125362**

### **BCAR3 (NM\_003567) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	BCAR3 (NM_003567) Human Untagged Clone
Tag:	Tag Free
Symbol:	BCAR3
Synonyms:	AND-34; MIG7; NSP2; SH2D3B
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)
Cell Selection:	None



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_003567, the custom clone sequence may differ by one or more nucleotides

```

ATGGCTGCAGGAAAATTGCAAGCCTCCAGAAACATGCCGGTGAATCACCAGTTCCTCCCTGGCCTCAT
CCATGGACCTTCTGAGCAGCAGGTCCCCTCTCGCTGAGCATCGCCAGATGCCTATCAAGATGTGTCTAT
ACATGGCACCTTCCACGGAAGAAAAAGTCTCTCCATAAGTCTGTGACTTCACTCAGTCACATG
GGCACCTCCCCACTCAAATCCCACGGCAGAACTCGCTGTGACCCAGGATGGCATCCAGGAGAGCC
CATGGCAGGACCGCACGGCGAAACCTTCACTTCAGGGATCCACATCTTCTGGACCCAACCTGTGGAATA
TGTGAAGTTCTCAAGGAGAGGCACATCATGGACAGGACCCCGAGAACTGAAGAAGGAGCTGGAGGAG
GAGCTGCTCTGAGCAGCAGGACCTGCGCAGCCATGCCTGGTACCACGGCCGCATCCCCGACAGGTGT
CTGAAAACCTTGTGCAGCGAGATGGTGACTTCTAGTTCGTGACTCTCTGTCCAGCCCTGGGAACCTTGT
CCTGACCTGTGAGTGAAGAACCTCGCTCAGCACTTCAAATCAACCGGACAGTTCTGCGACTCAGCGAG
GCCTACAGCCGCTGCAGTACCAGTTCGAGATGGAGAGCTTCGACTCCATCCCCGGCCTGGTGCCTGCT
ACGTGGGCAACCGCCGGCCATCTCCAGCAGAGTGGCGCCATCATCTTCCAGCCATCAACAGGACGGT
GCCTCTGCGGTGCTGGAGGAGCATTATGGCACCTCCCCAGGCCAGGCCGGGAGGGCAGCTCACCAGG
GGAAGGCCGATGTGGCAAGAGGCTGAGCCTCACCATGGGTGGCGTCCAGGCCGAGAGCAGAATTTGC
CCAGGGGAAACCTCCTCAGAAACAAAGAAAAGAGTGGTAGCCAGCCCGCCTGCCTGGATCAGATGCAGGA
CAGAAGAGCCTTGTCCCTCAAAGCCCACAGTCAGAGAGCTACCTGCCGATTGGCTGCAAGCTGCCACCT
CAGTCTCGGGTGTGGACACAAGCCCCTGCCAAACTCACCTGTGTTTCCAGGACGGGAAGCGAGCCTGCC
TGAGCCCAGCAGTGGTTCGGAGGGTCTCCTCAGACGCCAGGGCTGGGGAGGGCTGAGGGGATCAGACAG
TCAACTGTGCCCTAAGCCCCGCCTAAGCCCTGCAAGGTGCCGTCTCTCAAGTTCCCTCGTCTCCCTCT
GCCTGGCTCAACTCAGAGGCCAACTACTGTGAAGTGAACCCAGCGTTTCCACAGGCTGCCGAGGGGAG
CAAAGCTACCCTCATGTGCCAGGGGAGCCACACAGAAGTCTCACAGCCAAGCAGAATGAGGCCCGCAGG
TCCCCGAACTCTGGCGTCAACTACTTGATCCTTGATGATGATGACAGGGAAAGACCTTGGGAACCTGCG
GCAGCTCAGATGGAGAAGGGGAGTGGGACAAGGGCGAGTTTGTGACGCCCTCCTGGAGACTGTCTCCT
CCTTCAGGCCAACGAGTTTGTGATCAAAGTTCTTCCCCCTGAGAATAAGCCCCTGGAAACAGCAATGTT
GAAACGTGCAAAAGAACTGTTACCAACAACGACCCCAAGGTCATCGCCAGCACGTAAGTACTGAGCATGGAC
TGCAGGGTTGCTAGGATACTTGGAGTCTCTGAAGAGATGAGGAGGAACATGGGGGTGAGCTCAGGCCTGG
AACTCATTACCTTGCCCTACGGACACCAGCTGCGCCTGGACATAATTGAAAGACACAACACAATGGCCAT
CGGCATTGCAGTGACATTCTGGGATGCACGGGCACTTTGGAGGACCGAGCGGCCACTCTGAGTAAGATC
ATCCAGGTGGCGGTGAACTGAAGGATTCATGGGGACCTCTATTCTTCTCAGCTCTCATGAAAGCCC
TGGAAATGCCACAGATCAAAAGTTAGAAAAGACGTGGACTGCTCTGCGGCACCAGTACACCCAACTGC
CATTCTCTATGAGAAACAGCTGAAGCCCTTCCAGAACTCCTGCATGAAGGCAGAGAGTCCACATGTGTT
CCCCAAACAATGTATCAGTCCCCTGCTGATGCCGCTTGTGACGTTAATGGAGCGCCAGGCTGTGACTT
TTGAAGGAACCGACATGTGGGAAAAAACGACCAGAGCTGTGAAATCATGCTGAACCATTTGGCAACAGC
GCGATTTCATGGCCGAGGCTGCAGACAGTACCGGATGAATGCTGAGAGGATCCTGGCAGGTTTTCAACCA
GATGAAGAAATGAATGAAATCTGCAAGACTGAATTTCAAATGCGATTGCTATGGGGCAGCAAAGGTGCAC
AAGTCAATCAGACAGAGAGATATGAGAAATCAACCAGATTTAACTGCCCTCTCGCGTAAATTTGAAAC
TCCTCTGTAAAGCAGGCAGAGCTTTGA
    
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_003567

**Insert Size:** 3000 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_003567.3</a></u> , <u><a href="#">NP_003558.1</a></u>
<b>RefSeq Size:</b>	3203 bp
<b>RefSeq ORF:</b>	2478 bp
<b>Locus ID:</b>	8412
<b>UniProt ID:</b>	<u><a href="#">O75815</a></u>
<b>Cytogenetics:</b>	1p22.1
<b>Domains:</b>	SH2, RasGEF
<b>Protein Families:</b>	Druggable Genome
<b>Gene Summary:</b>	<p>Breast tumors are initially dependent on estrogens for growth and progression and can be inhibited by anti-estrogens such as tamoxifen. However, breast cancers progress to become anti-estrogen resistant. Breast cancer anti-estrogen resistance gene 3 was identified in the search for genes involved in the development of estrogen resistance. The gene encodes a component of intracellular signal transduction that causes estrogen-independent proliferation in human breast cancer cells. The protein contains a putative src homology 2 (SH2) domain, a hall mark of cellular tyrosine kinase signaling molecules, and is partly homologous to the cell division cycle protein CDC48. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2012]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein (isoform 1).</p>