

Product datasheet for SC125702

TRA (BC063385) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRA (BC063385) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRA
Synonyms:	TRA
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC063385 edited

```

CTGGGACTGAGGTGTTCTAAGCCTGGAAGACACAAAGCTGTCTCATCCCAACCTGCCG
TGCCCCCTCCTTGAGGGTTAGTGTAAGGCTCTGAAGGACTGTGTGAATTATGGTGGTGCTA
CAAACAAGCTCATCTTTGGAAGTGGCACTCTGCTTGCTGTCCAGCCAAGTACATATCCAG
AACCTGACCCTGCCGTGTACCACTGAGAGACTCTAAATCCAGTGACAAGTCTGTCTGC
CTATTCACCGATTTTGATTCTCAAAACAAATGTGTCAAAAGTAAGGATTCTGATGTGTAT
ATCACAGACAAAACCTGTCTAGACATGAGGTCTATGGACTTCAAGAGCAACAGTGTGTG
GCCTGGAGCAACAAATCTGACTTTGCATGTGCAACGCCTTCAACAACAGCATTATTCCA
GAAGACACCTTCTCCCCAGCCAGAAAGTTCTGTGATGTCAAGCTGGTCGAGAAAAGC
TTTGAAACAGATACGAACCTAACTTTCAAAACCTGTCAGTGATTGGGTTCCGAATCCTC
CTCCTGAAAGTGGCCGGGTTTAATCTGCTCATGACGCTGCGGCTGTGGTCCAGCTGAGAT
CTGCAAGATTGTAAGACAGCCTGTGCTCCCTCGCTCCTTCTCTGCATTGCCCTCTTCT
CCCTCTCCAAACAGAGGGAACCTCTCTACCCCCAAGGAGGTGAAAGCTGTACCACCTCT
GTGCCCCCCCCGGCAATGCCACCAACTGGATCCTACCCGAATTTATGATTAAGATTGCTGA
AGAGCTGCCAAACACTGCTGCCACCCCTCTGTTCCCTTATTGCTGCTTGTCACTGCCTG
ACATTCACGGCAGAGGCAAGGCTGCTGCAGCCTCCCTGGCTGTGCACATTCCCTCCTGC
TCCCCAGAGACTGCCTCCGCCATCCCACAGATGATGGATCTTCAGTGGGTTCTCTGGGC
TCTAGGTCCTGGAGAATGTTGTGAGGGGTTATTTTTTTTAATAGTGTTCAATAAGAAAT
ACATAGTATTCTTCTTCAAGACGTGGGGGAAATTATCTCATTATCGAGGCCCTGCTA
TGCTGTGTATCTGGCGTGTGATGTCTGCTGCCGATGCCTTCATTAATAATGATTTGG
AAGAGCAAAAAAAAAAAAAAAAAAAAAAAAAA
  
```

Restriction Sites:	NotI-NotI
ACCN:	BC063385



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC063385.1 , AAH63385.1
RefSeq Size:	1172 bp
Locus ID:	6955
Cytogenetics:	14q11.2