

Product datasheet for **SC125883**

SEC61G (BC009480) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SEC61G (BC009480) Human Untagged Clone
Tag:	Tag Free
Symbol:	SEC61G
Synonyms:	SSS1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC009480 edited CTGGCATTTTAGGTGTCGGTTGGGTAGGCAGTCATGGATCAGGTAATGCAGTTTGTGAG CCAAGTCGGCAGTTTGTAAAGGACTCCATTCGGCTGGTTAAAAGATGCACTAAACCTGAT AGAAAAGAATTCCAGAAGATTGCCATGGCAACAGCAATAGGATTTGCTATAATGGGATTC ATTGGCTCTTTGTGAAATTGATCCATATTCCTATTAATAACATCATTGTTGGTGGCTGA ATACATTTTGAAGAGAGTTTTTCATCTTAGAGATTGGTGAACAAGTGTGAGGGTGTGAG AAACTCACAGAATACAAATTTGCCTGTATGTTTTGTGGGTTTTTTTTTTCCTTTCAAGAT GTTTTCTATTTCTAAATTAAGTAATTTCAAAGTAAAAAAAAAAAAAAAAAAAAAAAAAAAA AAAAAAAAAAAAAAAAAAAA
5' Read Nucleotide Sequence:	>OriGene 5' read for BC009480 unedited NGAAGTTCGGATTTTGTAAACGACTTACTATAGGCGGACCGCAATTCGCACGAGGCTG GCATTTTAGGTGTCGGTTGGGTAGGCAGTCATGGATCAGGTAATGCAGTTTGTGAGCCA AGTCGGCAGTTTGTAAAGGACTCCATTCGGCTGGTTAAAAGATGCACTAAACCTGATAGA AAAGAATTCCAGAAGATTGCCATGGCAACAGCAATAGGATTTGCTATAATGGGATTCATT GGCTTCTTTGTGAAATTGATCCATATTCCTATTAATAACATCATTGTTGGTGGCTGAATA CATTTTGAAGAGAGTTTTTCATCTTAGAGATTGGTGAACAAGTGTGAGGGTGTGAGAAA CTCACAGAATACAAATTTGCCTGTATGTTTTGTGGGTTTTTTTTTTCCTTTCAAGATGTT TTCTATTTCTAAATTAAGTAATTTCAAAGTNNANAAAANNAANAAAAAAAAAAAAAAAA AAAAAAAAAAAAACCTCGACTCTAGATTGCGGCCGCGGTCATAGCTGTTTCCTGAACAAATC CCGGGTGGCATCCCTGTGACCCCTCCCCAGTGCCTCTCCTGGCCCTGGAAGTTGCCACTC CAGTGCCACAGCCTTGTCTAATAAAATTAAGTTGCATCATTTTGTCTGACTAGGTGT CCTTCTATAATATTATGGGGTGGAGGGGGTGGTATTGAACAAGGGGCAATTGGGAAAAA CAAACCGAGGGCCCTGGGGTCTATGGGAACCAGCTGGGGTGCAGTGGCCCAATCTGG GTTAATGGAATCTCCGCCTCCGGTTCAAGCGATTTCTGCCCCACCTCCCGATTTTGG GTTCCCGCCGCTAAACGGCTACCAATTTTTGTTTTTGGTAAAAAGGGTTTCCA
Restriction Sites:	Please inquire



[View online »](#)

ACCN:	BC009480
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC009480.2 , AAH09480.1
RefSeq Size:	436 bp
RefSeq ORF:	204 bp
Locus ID:	23480
Cytogenetics:	7p11.2
Protein Families:	Transmembrane
Protein Pathways:	Vibrio cholerae infection
Gene Summary:	The Sec61 complex is the central component of the protein translocation apparatus of the endoplasmic reticulum (ER) membrane. Oligomers of the Sec61 complex form a transmembrane channel where proteins are translocated across and integrated into the ER membrane. This complex consists of three membrane proteins- alpha, beta, and gamma. This gene encodes the gamma-subunit protein. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Jul 2008]