

Product datasheet for SC125930

RFPL2 (BC051910) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RFPL2 (BC051910) Human Untagged Clone
Tag:	Tag Free
Symbol:	RFPL2
Synonyms:	ret finger protein-like 2; RNF79
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC051910 edited

```

GCCAGCCCATCTGCTCCATGCCAGGCTTTTCTGGCTTTGTCTGCCCGGATGCCAGTGC
CTTGAGAAAGGAGCTACATGCCACGGGCTGGAGGCTAGAGCCTGTGGCACCATTGTTGCC
TCCGATGGAGGTGGCTGAATTAGGCTTCCAGAGACTGCAGTGTCCCAATCCAGGATCTG
TCTATGTGCTGTATTGTGTGGCACTGGGACTTTGCAGACATGATGGTGATAAGGAGCCT
GAGTGGACATGGCTGCACTCTTCCAAGAAGCAAGCAGCTGTCCCGTCTGCTCAGACTATC
TGGAAAAACCAATGTCCCTGGAGTGTGGATGCGCCGTCTGCCTCAAGTGCATTAATTCAC
TGCAGAAGGAGCCCCATGGGGAGGATCTACTTTGCTGTTGCTCTTCCATGGTCTCTCGGA
AGAACAAAATCAGGCGCAATCGGCAGCTAGAGAGGCTGGCTTCCACATCAAGGAACTGG
AGCCCAAGCTGAAGAAGATTCTGCAGATGAACCAAGGATGCGGAAGTTCAAGTGGATA
TGACCTTGGATGCCAACACAGCCAACAACCTTCTCCTCATTTCTGACGACCTCAGGAGCG
TCCGAAGTGGGCGCATCAGACAGAATCGGCAAGACCTTGCCGAGAGATTTGACGTGTCCG
TTTGCATCCTGGGCTCCCTCGCTTTACCTGTGGCCGCCACTGTGGGAGGTGGACGTGG
GAACAAGCACAGAATGGGACCTGGGAGTCTGCAGAGAATCTGTTACCGCAAAGGGAGGA
TCCAGCTGACCACAGAGCTTGGATTCTGGACTGTGAGTTTGAGGGATGGAGGCCGCTCT
CTGCCACCACGGTGCCGCTGACTTTCCTTCTCGTAGACCGCAAGTTACAGCGAGTGGGA
TTTTTCTGGATATGGGCATGCAGAACGTTTCTTTTTTGATGCTGAAAGTGGTCCCATG
TCTATACATTACAGGAGCGTATCTGCTGAGGAGCCATTGCGCCATTTTGGCTCCTTCAG
TTCCACCTAATGGTGATCAAGGTGTCTTGAGCATCTGTCTTTGATGAACTCAGGCACTA
CTGATGCTCCAGTCCGTCTGGGAGGCCAAATAGCCCTCACTCCAAAAAACAAAAAA
CAGGGTAAGAAAATTACTTGGGTGGGTAGACTTAGGAACGCTCTACTTCGTAAAAATCATT
ATACAAAGTCACGGGAGAAAAATATGGGACATTTCTTGATTGTACTTAATCTAATTTGAT
TAGATTATAGAGTCCTAAGTATTAATTATTGCCACCATCAAACCTATTGAGTCCTATGGT
TCACATCTTGTTCCTATAGAAATGTCCTGTATTCTGGGATCAATTTCAAATGCTTTAC
TTTTTTATTTCTGCAAGTTCAAATTAATGTATTATAGAAGTTATGAGTTAAATAAAGATT
GGAATATCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAA
  
```



[View online »](#)

Restriction Sites:	NotI-NotI
ACCN:	BC051910
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC051910.1</u> , <u>AAH51910.1</u>
RefSeq Size:	1515 bp
Locus ID:	10739
Cytogenetics:	22q12.3
Protein Families:	Druggable Genome