

## Product datasheet for **SC125994**

### UQCRH (BC001426) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	UQCRH (BC001426) Human Untagged Clone
Tag:	Tag Free
Symbol:	UQCRH
Synonyms:	MGC111572; OTTHUMP0000009481; QCR6; ubiquinol-cytochrome c reductase hinge protein
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL5</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for BC001426 edited  
 CTGGGTTAGGTGCCGCTGTTGCTGCTCGTGTGAATCTAGAACCGTAGCCAGACATGGGA  
 CTGGAGGACGAGCAAAAGATGCTTACCGAATCCGGAGATCCTGAGGAGGAGGAAGAGGAA  
 GAGGAGGAATTAGTGGATCCCCTAACAAACAGTGAGAGAGCAATGCGAGCAGTTGGAGAAA  
 TGTGTAAGGCCCGGAGCGGCTAGAGCTCTGTGATGAGCGTGTATCCTCTCGATCACAT  
 ACAGAAGAGGATTGCACGGAGGAGCTCTTTGACTTCTTGCATGCGAGGGACCATTGCGTG  
 GCCCACAAACTCTTTAACAACTTAAAATAAATGTGTGGACTTAATTCACCCCAGTCTTCA  
 TCATCTGGGCATCAGAATATTTCTTATGGTTTTGGATGTACCATTTGTTTCTATTTGT  
 GTAAGTTACATGAACCTCATGGGTTTGGCTTAGGCTGGTAGCTTCTATGTAAT  
 TCGCAATGATTCCATCTAAATAAAAGTTCTATGATCTGCAAAAAAAAAAAAAAAAAAAAA



[View online »](#)



**Gene Summary:**

This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This protein may mediate formation of the complex between cytochromes c and c1.[UniProtKB/Swiss-Prot Function]