

Product datasheet for **SC126022**

ZNF695 (BC041082) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF695 (BC041082) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF695
Synonyms:	SBZF3
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC041082 edited

```

CTGGGTTCCCGTCTGCTGCCTCTCGGAGAGTCCCGGGTGACTGCCGCAGGCTCCATCGCC
CTGTGGCCTGCAGGTATTGCGAGATTTATAGGGAGGACGCTGGGACCCCAAAAGCTGGG
AAATGGGACTATTGGCATTACAGGATGTGGCTCTAGAATTCTCTCCAGAGGAGTGGAAT
GCCTGGACCCAGCTCAGCGGAGTTTGTATAGGGATGTGATGTTAGAGAACTACAGAAACC
TGATCTCCCTTGGTGAGGATAGCTTCAATATGCAATTCCTATTTACAGTCTTGCTATGT
CTAAGCCAGAAGTATCATCTGTCTGGAGGCAAGGAAAGAGCCCTGGAACGTGAACACAG
AGAAGACAGCCAAACACTCAGCTTTGTCTTCTTACTTACTGAAGACATTTGCCAGAGC
AGGGCCTGCAAGTTTCATTCCAAAAAGTGATGCTGAGAAGATATGAAAGATGTTGCTTG
AGAAATTACGCTTAAGGAATGACTGGGAAATTCATGTGAAGATGTGCTTGCTTCCCTT
TGCTTCTGCCATGATTCTAAGTTTCCTGAGGCCTCCCCAGAAGCAGAAGCATGTAAAGC
CCACAGAACCGGTCCAGTCGTCGACCATCGAGTGCTTTGAAGACCTGCGGAGGTCCG
TCGCTGCTTGAACACAGTGCCTGGACTCCCTGCTCCGTGAGGCGAATTTACACCGGGTCC
AGCTGGGTGAGCCCGAGTTCTGAATTAACATGCCTCCCAAAAAAAAAAAAAAAAAAAAA
  
```


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC041082 unedited <pre> GTTTAGATTTGTATACGACTTATATAGGGCGGCCGCAATTTCGCACGAGGCTGGGTTCCC GTCTGCTGCCTCTCGGAGAGTCCCGGTGACTGCCGCAGGCTCCATCGCCCTGTGGCCTG CAGGTATTGCGAGATTTATAGGGAGGACGCTGGGACCCCCAAAAGCTGGGAAATGGGACT ATTGGCATTACAGGATGTGGCTCTAGAATTCTCTCCAGAGGAGTGGGAATGCCTGGACCC AGCTCAGCGGAGTTTGTATAGGGATGTGATGTTAGAGAACTACAGAAACCTGATCTCCCT TGGTGAGGATAGCTTCAATATGCAATTCCTATTTACAGTCTTGCTATGTCTAAGCCAGA ACTGATCATCTGTCTGGAGGCAAGGAAAGAGCCCTGGAACGTGAACACAGAGAAGACAGC CAAACACTCAGCTTTGTCTTCTTATCTTACTGAAGACATTTTGCCAGAGCAGGGCTGCA AGTTTCATTCCAAAAAGTGATGCTGAGAAGATATGAAAGATGTTGTCTTGAGAAATTACG CTTAAGGAATGACTGGGAAATCCATGTGAAGATGTGCTTGCTTCCCTTTGCCTTCTGC CATGATTCTAAGTTTCTGAGGCCTCCCAGAAGCAGAAGCATGTAAAGCCACAGAACC GGTCCAGTCGTCGACCATCGAGCTGCTTTGAAGACCCTGCNGAGGTTCCGTCGCTGCTTT GACACAGTGCCTGGACTCCCTGCTNCGTGAGGGCGATTTACACCGGGTCCGCTGGGGTC AGCCCCGAGTTCTGAAATAAAACATGCCTNCCACAAAAAATACTCGACTCTA GATGGCGCCGCCGTCATAGCTGTTTCTGAAAGATCCCGTTGCATCCCTGTGCCCTCC CCATGCTCTCCTGCCCTGGAAGTGCCCTCC </pre>
Restriction Sites:	Please inquire
ACCN:	BC041082
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC041082.1 , AAH41082.2
RefSeq Size:	780 bp
Locus ID:	57116
Cytogenetics:	1q44
Protein Families:	Transcription Factors
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]