

Product datasheet for SC126200

RPL14 (BC029036) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPL14 (BC029036) Human Untagged Clone
Tag: Tag Free
Symbol: RPL14
Synonyms: CAG-ISL-7; CTG-B33; hRPL14; L14; RL14
Mammalian Cell Selection: None
Vector: pCMV6-XL5
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for BC029036 edited
CGCCGCCAACATGGTGAGTCTTACTGTTGCGGGCTCCGGGCGTCGACCATGCCGCTCG
ACCTCCACCTCCGCTGGGAAGCTGAGGCGCCGAACGGCTCCAGAGGGTCCCGGAAGCG
CATGGTGTTCAGGCGCTTCGTGGAGTTGGCCGGTGGCCTATGTCTCCTTTGGACCTCA
TGCCGGAAAATTGGTCGCGATTGTAGATGTTATTGATCAGAACAGGGCTTTGGTCGATGG
ACCTTGCACTCAAGTGAGGAGACAGGCCATGCCTTTCAAGTGCATGCAGCTCACTGATTT
CATCCTCAAGTTCCGCACAGTGCCACCAGAAGTATGTCGACAAGCCTGCCAGAAGGC
AGACATCAATACAAAATGGGCAGCCACACGATGGGCCAAGAAGATTGAAGCCAGAGAAAG
GAAAGCCAAGATGACAGATTTTGTATCGTTTTAAAGTTATGAAGGCAAAGAAAATGAGGAA
CAGAATAATCAAGAATGAAGTTAAGAAGCTTCAAAGGCAGCTCTCCTGAAAGCTTCTCC
CAAAAAGCACCTGGTACTAAGGGTACTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTGCTG
TGCTAAAGTTCCAGCAAAAAGATCACCGCCGCGAGTAAAAGGCTCCAGCCAGAAAGGT
TCCTGCCAGAAAAGCCACAGGCCAGAAAGCAGCGCCTGCTCCAAAAGCTCAGAAGGGTCA
AAAAGTCCAGCCAGAAAGCACCTGCTCAAAGGCATCTGGCAAGAAAGCATAAGTGCC
AATCATAAAAAGTATAAAGGTTCTTTTTGACCTGTTCAAAAAAAAAAAAAAAAAAAAAA
AAAAAAAAAAAAAAAAAAAAA



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC029036 unedited GGCCGTTCAAATTTTGTATACGACTCATATAGGCGGCCGCGNAATCANATCTGGTACCGG TCCGGAATCCCGGGATCGCCGCCAACTGGTGAGTCTTACTGTTGCGGGCTCCGGGGCCG TCGACCATGCCGCTCGACCTCCACCTCCGCTGGGAAGCTGAGGCGCCGAACGGCTCCCAG AGGGTCCCAGGGAAGCGCATGGTGTTCAGGCGCTTCGTGGAGGTTGGCCGGGTGGCCTATG TCTCCTTTGGACCTCATGCCGAAAATTGGTCGCGATTGTAGATGTTATTGATCAGAACA GGGCTTTGGTCGATGGACCTTGCACTCAAGTGAGGAGACAGGCCATGCCTTTCAAGTGCA TGCAGCTCACTGATTTTCATCCTCAAGTTTCCGCACAGTGCCACCAGAAGTATGTCCGAC AAGCCTGGCAGAAGGCAGACATCAATACAAAATGGGCAGCCACACGATGGGCCAAGAAGA TTGAAGCCAGAGAAAGGAAAGCCAAGATGACAGATTTTGTATCGTTTTAAAGTTATGAAGG CAAAGAAAATGAGGAACAGAATAATCAAGAATGAAGTTAAGAAGCTTCAAAGGCCAGCTC TCCTGAAAGCTTCTCCCAAAAAGCACCTGGTACTAAGGGTACTGCTGCTGCTGCTGCTG CTGCTGCTGCTGCTGCTGCTAAAGTTCCAGCAAAAAGATCACCGCCGCGAGTAAAAAGG CTCCAGCCAGAAAGTTCTGCCCCAGAAGCCACAGGCCAGAAAGCAGCGCCTGCNTCAA AAGCTCAGAAGGGTCAAAAAGCTCCAGCCAGAAAGCACCTGCTCAAAGCATCTGGCCA GAAAGCATAGTGGCATCATAAAAGTAATAAAGGTCCTTTTGACCTGTTTCAAAAAAAAAA AAAAAA
Restriction Sites:	Please inquire
ACCN:	BC029036
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC029036.1 , AAH29036.1
RefSeq Size:	858 bp
Locus ID:	9045
Cytogenetics:	3p22.1
Protein Pathways:	Ribosome

Gene Summary:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L14E family of ribosomal proteins. It contains a basic region-leucine zipper (bZIP)-like domain. The protein is located in the cytoplasm. This gene contains a trinucleotide (GCT) repeat tract whose length is highly polymorphic; these triplet repeats result in a stretch of alanine residues in the encoded protein. Transcript variants utilizing alternative polyA signals and alternative 5'-terminal exons exist but all encode the same protein. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]