

Product datasheet for **SC126298**

CSAG1 (NM_153478) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CSAG1 (NM_153478) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSAG1
Synonyms:	CSAGE; CT24.1
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_153478 edited CTGGATTCTTCTGGATGGGGATCCAGATGGAGGTGGAGGTTGATTTGGGAAGCAGAGC ACAGCAGCCCAAATTTGCTTGTAAATGTCGGCGACTACAGCCTGCTGGCCTGCCTTCACTG TCCTGGGGGAAGCTCGGGGAGACCAGGTGGACTGGAGTAGACTGTACAGAGACACTGGTC TGGTGAAGATGTCCAGGAAACCACGAGCCTCCAGCCCATTTTCCAACAACCAACCATCAA CACCAAAGAGGTTCCCAAGACAACCCAGAAGGGAAGGGACCCGTCAGGAAGTTCAG GAACAAAAGGCTCTCCCTAAAAGACCACCGCTTCAAAAAACCTGAGGAATGGAGTGGGC CAACACTATCCAGCCACTCTGACCAGCCGAACGAGGAAGTCAATCAAAATGCGCCATAGC AGGACCACAAGGGCAAGGAGACCACCGCTTCTCCAGTGCTTCTTGGGCAGCCAGTAAT TCCCAGGCAAGGCCAGAGACTTCAAGTCTATCTGAAAAGTCTCCAGAAGTCTAACCCAG ATAAATAGCCAACAGGGTGTAGGTACGTTTTACACCCAAAGGGTATGCCCCATGGTGA TGGAATAAAATGAACATGTTGTAAAAAAAAAAAAAAAAAAAAA



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for NM_153478 unedited <pre> CGGTTCATATTTGTATACGACTCATATAGGCGGCCGCGTATTCAGATCTGGTACCGGTCC GGAATTCGCGGATATCGTCGACCCACGCGTCCGCCACGCGTCCGCTGGATTCTTCCTG GATGGGGATCCAGATGGAGGTGGAGGGTTGATTGGGAAGCAGAGCACAGCAGCCAAAT TTGCTTGAATGTCGGCGACTACAGCCTGCTGGCCTGCCTTCACTGTCTGGGGGAAGCT CGGGGAGACCAGGTGGACTGGAGTAGACTGTACAGAGACACTGGTCTGGTGAAGATGTCC AGGAAACCACGAGCCTCCAGCCATTTTCCAACAACCAACCATCAACACCAAGAGGTTT CCAAGACAACCCAGAAGGAAAAAGGACCGTCAAGGAAGTTCCAGGAACAAAAGGCTCT CCCTAAAAGACCACCGCTTCAAAAAAACCTGAGGAATGGAGTGGGCAACACTATCCAGC CACTCTGACCAGCCGAACGAGGAAGTCAATCAAAATGCGCCATAGCAGGACCACAAGGGC AAGGAGACCACCGCCTTCTCCAGTGCTTCTTGGGCAGCCAGTAATCCAGGCAAGGCC AGAGACTTCAAGTCTATCTGAAAAGTCTCCAGAAGTCTAACCCAGATAAATAGCCAACA GGGTGTAGAGTACGTTTTACACCCANAGGGTATGCCCCATGGTGATGGAATAAAATGA ACATGTTGTAAAAAAAAAAAAAAAAAGGGCGGCCGCGGTATAGCTGTTTCTGAAC AGATCCCGGGTGGCATCCCTGTGACCCCTCCAGTGCCTCTCCTGGCCCTGGAAGTGCC ACTCCAGTGCCCAACAGCCTTGTCTAATAAAATTAAAGTGCATATTTGTCTGACTAG GTGTCCTTCTATAATATATGGGGTGAGGGGGTGGGAAA </pre>
Restriction Sites:	Please inquire
ACCN:	NM_153478
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_153478.1</u> , <u>NP_705611.1</u>
RefSeq Size:	885 bp
RefSeq ORF:	237 bp
Locus ID:	158511
UniProt ID:	<u>Q6PB30</u>
Cytogenetics:	Xq28

Gene Summary:

This gene encodes a member of a family of tumor antigens. The protein is expressed in chondrosarcomas, but may also be expressed in normal tissues such as testis. Alternative splicing of this gene results in two transcript variants encoding the same protein. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (a) represents the longer transcript. Variants a and c both encode the same protein.