

Product datasheet for **SC126300**

C16orf74 (BC009078) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C16orf74 (BC009078) Human Untagged Clone
Tag:	Tag Free
Symbol:	C16orf74
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC009078 edited GCAGCCAGATCTGCTGGGACACCTTTCCCAAGGAAGAGCCCGTTGCACTGGGCTTTGAAG GATAAGCAGGAGCTTGTACTCAGGCAGAGGAAGAAAGAGCATCCCAGGCGGGGGAGCA GCATATGCAAAGGCACGAAGGGGCCCCAGGAGCCTAGGGAGTCTGGGGAAGTGTGAGCAC TTTGGAGAGTGGAGGCTGGAGCGCTGTGGAGAGTGGGGGCTGGTGGCCGGAATGAAGCT GCAGCTGGCTGGGCCACATGGTAAAGGCTGACAACTGGACCCAGAGGCCAACTAGCCTAT GATCAGCATTTCCCAAAATCTGTTTCCGACTCATGGTTCTGTGAGATGTGACAAGGGCT CCTTTTTTCATTCTGAGACGCCGGTTTTTCATCTGTGATGCGGGGACAGCTGCGCTCCTTG CTGCGAGGCGTCAGGACCCAGGTGATAGTGAAGGGAGGGTGGCGCCCGCGGTTCCCGGCG GCCACTGATGCCTGTCTCTCTGTGCTGTGTACGTGCGTGTGTGCTCCACGCCTGGCTTCT CAGGCTTTCAAATGTGTGTCAGCAGCAGCAGCAGCAGCCAGCAGGCCCCCGTCTGA ACGACAAGCACCTGGACGTGCCCGACATCATCATCACGCCCCCACCACCGGCATGA TGCTGCCGAGGGACTTGGGGAGCACAGTCTGGCTGGATGAGACAGGGTCGTGCCAGATG ATGGAGAAATCGACCCAGAAGCCTGAGGAGGTGTCTGGGTTTGGCTGGCTGGCTCCTGC TCCAGCGGCCCGGCTTCAGGTGTCCGGGGCGTGGCTGCCTGGAGCAGGTGTGCTGAATA CCCTGGATGGGAAGTGAAGCAACCCGGGCTCCGCTCAGAGAGACGTGGCAGGACCAGCG AGGAATCCAGCCTGTCCACTTCCAGAACAGTGTTCAGGCCCCGCTGAGTGGACCGGA CCTCTGACACCTCAGGTTCTTGCTGACTCCGGCCTGGTAAAGGGAGCGCCATGGTCTCT GGCTGTTGGGGTCCCAGGAGAGGCTCTCTTCTGGACAAACACACCTCCAGCCCCAG GGCTGTGCAAACACATGCCCTCCCATAGCACCAACAAGAACTTCTTGCAGGTGGAGTG GCTGTTTTTTATAAGTTGTTTTACAGATACGGAACAGTCCAAAATGGGATTTATAATTT CTTTTTTGCAATTAAATAAAGATCCTCTGTAAACAAAAAAAAAAAAAAAAAAAAAAAAA A
Restriction Sites:	Please inquire
ACCN:	BC009078



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>BC009078.1</u>
RefSeq Size:	1261 bp
Locus ID:	404550
Cytogenetics:	16q24.1