

Product datasheet for **SC126314**

C3orf14 (BC017772) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C3orf14 (BC017772) Human Untagged Clone
Tag:	Tag Free
Symbol:	C3orf14
Synonyms:	chromosome 3 open reading frame 14; HT021
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC017772 edited

```

GGCAAGCACTTTAACCTTTTAAGCCCAACCAGATGAGTTGCCTGCAGTTTTGGAGGCCTT
CAGAGCATTTCACTAGACCTCTGTCTGTGTCGGTCCAGTGTCTTTAGCCAAGCTTTGATT
AAAGATGACTTCCTTGTTTGCTCAAGAAATTCGCCTTTCTAAAAGACATGAAGAAATAGT
ATCACAAGATTAATGTTACTTCAACAAATGGAGAATAAATTGGGTGATCAACACACAGA
AAAGGCATCTCAACTCCAACTGTTGAGACTGCTTTTAAAGGAACCTTAGTCTTTTAA
GGATATAGAAGCAGCAGAAAAGTCACTACAGACCAGGATTCACCACTTCCACGGCCTGA
GGTGGTTTCTCTTGAGACTCGTTACTGGGCATCAGTAGAAGAATATATCCCAAATGGGA
ACAGTTTCTTTTAGGAAGAGCACCATATCCTTTTGCTGTTGAAAATCAAATGAAGCAGA
AAATACCATTCAAATGAGGCACAGCGATAACTTCTTCACATGCTATTTCAAAAAGCCTG
TTTAATAAAGCTGAATGTTAAGGTGTATGTAGGTATTGCAGGAACCTTAGGAATTAAT
ATGTTTCATATTCTCGATTATCTCCTAAGTGACAGTGAAGATATGAGAATTTACTGGCAA
GTCACATGTTATCACCTACTACTATTTCAAGGTCATGAATTTGCTTTCTACCAAAACCCAT
AGATGTGTTAAACACGAATATTAAGGTGGACTCTTTATTAATTAACCAAAACCCAT
AAAAAAAAAAAAAAAAAAAAACCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA

```



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC017772 unedited <pre> CCGGGGGTCAGATATTTGTATACCACTCACTATAGGCGGCGCGGATTGCGCATTACGG GCCGGGGGCTTCCCGGAGGGCTCAGAGGCGGGCTCGGGCAAGCACTTTAACCTTTAAG CCCAACCAGATGAGTTGCCTGCAGTTTTGGAGGCCTTCAGAGCATTCTACTAGACCTCTG TCTGTGTCGGTCCAGTGTCTTTAGCCAAGCTTTGATTAAGATGACTTCCTTGTTTGCTC AAGAAATTCGCCTTTCTAAAAGACATGAAGAAATAGTATCACAAGATTAATGTTACTTC AACAAATGGAGAATAAATTGGGTGATCAACACACAGAAAAGGCGTCTCAACTCCAACTG GTGAGACTGCTTTTAAAGGAACCTTAGTCTTTTAAAGGAGATAGAAGCAGCAGAAAAGT CACTACAGACCAGGATTACCCACTTCCACGGCCTGAGGTGGTTTCTCTTGAGACTCGTT ACTGGGCATCAGTAGAAGAATATATTCCCAATGGGAACAGTTTCTTTAGGAAGAGCAC CATATCCTTTTGTGTTGAAAATCAAAATGAAGCAGAGAATACCATTCAAAATGAGGCAC AGCGATAACTTCTTCACATGGCTATTTCAAAAAGCCTGTGAATAAAAGCTGAATGTAA GGTGTATGTAGGTTATTGCAGGGACTTTTAGAATTAATATGTTTATCTTCGATTAT CTCCCAAGTGACAGGGGAAATATGAGAATTTACTGGCGGGCAATGGTTATCCCTAATAC TATTTCCAGGTGAGGAATTTGCTTTCACCAACCCTAAGGGTGGTAAACCCGATATTA AGAGGGGGCGCTCTAATTAAGTTGGGACCGAGAAGAGGGATTGAACGCGAT </pre>
Restriction Sites:	Please inquire
ACCN:	BC017772
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC017772.2 , AAH17772.1
RefSeq Size:	832 bp
Locus ID:	57415
Cytogenetics:	3p14.2