

Product datasheet for **SC126334**

FLJ40125 (BC062452) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FLJ40125 (BC062452) Human Untagged Clone
Tag:	Tag Free
Symbol:	FLJ40125
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for BC062452 edited <pre> GGGGAGAGTTAGTGGAAGGGACTTAAGAGAAAGGGCGTGGTCTTTTGGAAAGGGGCGATT CTGGGTCGTAGGAGCCGGCCGGAAGTATTGATTAGTGCTGGAAGGTGGAAGAGCAGGGC AGCCTGGACAACATGACCTGCATCCTGGTCTGCTTCCCTGGGGCCCTAGGCCTTCTGAG GAGGCGATCAGGAGGGAGCTAGCACTGGACGCAGCCCTGGGCTGCAGAATCGCTGAACTG TGTGCCTCTGCTCAGAAGCCCCCAGCCTGAACACAGTTTTTCAGGACTCTGGCCTCAGAG GACATCCCAGATTTACCTCCTGGGGGAGGGCTGGACTGCAAGGCCACTGTCATTGCTGAA GTTTATTCTCAGATCTGCCAGGTCTCAGAAGAGTGCGGAGAGAAGGGGCAGGATGGGGCT GGGAAGTCCAACCCACGCATTTGGGCTCAGCCTTGACATGGAGGCCTGACAGCTGTTG TCCTTTGGGGATCCTTTGCTTCTCTGGGGCCTCAACAGAACTAAAGAAGAAAACCGACCC TTTCCCCAACTACATGTACCAGCGGAAGGAAGGAAGGCCAATGTAGGAACCCAAAATGCT TATTTCTTCTTCTTACTTCCCTCTCACAGAAAAGTCTTACGAATGGGGAAATTCACC AACATCCAGACCAAAAAGAAAAAGCCCCAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA </pre>



[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC062452 unedited <pre> AGGTTCATATTTGTATACGACTCATATAGGCGGCCGCGTATTCGGCCATTACGGCCGGG AGAGTTAGTGGAAGGGACTTAAGAGAAGGGCGTGGTCTTTTGGAAAGGGGCGATTCTGGG TCGTAGGAGCCGGGCCGGAAGTATTGATTAGTGCTGGAAGGTGGAAGAGCAGGGCAGCCT GGACAACATGACCTGCATCCTGGTCTGCTTCCCTGGGGCCCTAGGCCTTCTGAGGAGGC GATCAGGAGGGAGCTAGCACTGGACGCAGCCCTGGGCTGCAGAATCGCTGAACTGTGTGC CTCTGCTCAGAAGCCCCCAGCCTGAACACAGTTTTTCAGGACTCTGGCCTCAGAGGACAT CCCAGATTACCTCCTGGGGGAGGGCTGGACTGCAAGGCCACTGTCATTGCTGAAGTTTA TTCTCAGATCTGCCAGGTCTCANAAGAGTGCGGAGAGAAGGGGCAGGATGGGGCTGGGAA GTCCAACCCACGCATTTGGGCTCAGCCTTGGACATGGAGGCCTGACAGCTGTTGTCTT TGGGGATCCTTTGCTTCTCTGGGCTCAACAGAACTAAAGAAGAAAACCGACCTTTCC CCAACTACATGTACCANCGGAAGGAAGGAAGGCCAATGTAGGAACCCAAAATGCTATTT CTTCTTCTTACTTCCCTCTCACAGAAAAGTCTTACGAATGGGGAAATCCACCAACAT CCAGACCAAAAAGGAAAAAGCCANNNNANANNNNNANNNNNAAAAAATGTTG GCCGCCTCGGCCCTCGCCTCTAGATTGCCGCCCGCGGTATAGCTGTTTCTGAACAGAT CCCGGGTGGCATCCCTGGGACCCCTTCCAGTGCCTCTCTTGCCCTGGGAGTGCCAC TCCGTGCCCCACCAACCTTTGCTCTAATAAAATAAAGTG </pre>
Restriction Sites:	NotI-NotI
ACCN:	BC062452
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC062452.1 , AAH62452.1
RefSeq Size:	719 bp
Locus ID:	147699
Cytogenetics:	19q13.32