

## Product datasheet for **SC127574**

### FCF1 (NM\_015962) Human Untagged Clone

#### Product data:

|                           |                                                                                                   |
|---------------------------|---------------------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                               |
| Product Name:             | FCF1 (NM_015962) Human Untagged Clone                                                             |
| Tag:                      | Tag Free                                                                                          |
| Symbol:                   | FCF1                                                                                              |
| Synonyms:                 | Bka; C14orf111; CGI-35; UTP24                                                                     |
| Mammalian Cell Selection: | None                                                                                              |
| Vector:                   | <u>pCMV6-XL4</u>                                                                                  |
| E. coli Selection:        | Ampicillin (100 ug/mL)                                                                            |
| Fully Sequenced ORF:      | >NCBI ORF sequence for NM_015962, the custom clone sequence may differ by one or more nucleotides |

ATGGGGAAGCAAAAGAAAACAAGGAAGTATGCGACCATGAAGCGAATGCTTAGTCTCAGAGATCAGAGGC  
TTAAAGAAAAGGATAGATTAAACCTAAAAAGAAAGAAAAGAAGGATCCCAGCGCATTAAAGGAAAAGAGA  
AGTTCCCCAACCCCTTCCTGCTTATTTTCCAATATAATACACAGCTGGGCCACCTTACCACATCCTC  
GTTGATACCAACTTTATCAACTTTCCATAAAAGCCAACTGGACTTAGTGCAGTCAATGATGGACTGTC  
TGTATGCCAAGTGTATCCCATGTATAACCGATTGTGTAATGGCTGAAATTGAGAAATTGGGCGAGAAGTA  
TCGAGTGGCTCTAAGGATTGCCAAGGATCCAAGATTTGAACGATTACCATGTACACACAAAGGAACCTAT  
GCAGATGACTGCTTAGTACAGAGAGTAACCTCAGCATAAGTGTTACATTGTGGCCACAGTTGACCGGGACC  
TGAAAAGAAGAATCCGTAAGATTCTGGAGTTCCTATCATGTACATTTCTAACCATAGATACAACATTGA  
ACGGATGCCAGATGATTATGGAGCCCCTCGATTCTAA



[View online »](#)

|                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>5' Read Nucleotide Sequence:</b> | <p>&gt;OriGene 5' read for NM_015962 unedited</p> <p>NAAGGCCAGAATTTTGTAAACGAACTCACTATAGGGCGGCCGCGATTCCGGCACCAGGGG<br/> AAAGCAAAAGAAAACAAGAAGTATGCGACCATGAAGCGAATGCTTAGTCTCAGAGATCAG<br/> AGGCTTAAAGAAAAGGATAGATTAACCTAAAAAGAAAAGAAAGGATCCCAGCGCA<br/> TTAAAGGAAAGAGAAGTTCCCAACACCCTTCTGCTTATTTTCCAATATAATACACAG<br/> CTGGGCCACCTTACCACATCCTCGTTGATACCAACTTTATCAACTTTTCCATAAAAGCC<br/> AACTGGACTTAGTGCACTCAATGATGGACTGTCTGTATGCCAAGTGTATCCCATGTATA<br/> ACCGATTGTGTAATGGCTGAAATTGAGAAATTGGGGCAGAAGTATCGAGTGGCTCTAAGG<br/> ATTGCCAAGGATCCAAGATTTGAACGATTACCATGTACACACAAAGGAACCTATGCAGAT<br/> GACTGCTTAGTACAGAGAGTAACCTCAGCATAAGTGTTACATTGTGGCCACAGTTGACCGG<br/> GACCTGAAAAGAAGAATCCGTAAGATTCTGGAGTTCTCTATCATGTACATTTCTAACCAT<br/> AGATACAACATTGAACGGATGCCAGATGATTATGGAGCCCTCGATTCTAATTCTTACAA<br/> GACACAGTTCTCTGCCTTTCTTCGACCAACTTTCTCTTGTGCCAGTTCATTACACANA<br/> ATGTAGCGGGATTNTTAAGGAATCAGAGAGACTGATGGNAGTCAGGNAGATATTTATTAT<br/> TTANGTGACCCAGCCAGTCAGATTAACATCCAAAGGACTGAACCCCTGACAGAGNTAAG<br/> NTACCTTTTAAGCATTTNNTGTGNGCCGCGNGGGGTNGNGGGGGATCTGTGCAGNGGAAG<br/> CATATTACAGAAGCAGAAC</p> |
| <b>3' Read Nucleotide Sequence:</b> | <p>&gt;OriGene 3' read for NM_015962 unedited</p> <p>TTTTTTTTTTTTTTTATTGTAAATTGACAATTTATAATTGTATAAATGTATGGGTACAAA<br/> GTGATATCATAAATTGTGGATACAATATGGAATAGTTAAATCAAGCTAGTTAACATATGC<br/> AATACTTCAAATACATTTCTGTGATGAGAACATCAATTTTGAAATGTATAATACCCATT<br/> ATTAAGTATAGGCACCACACTGTACAACAGAACTCAAAAAAGCATAAACATATTCTTC<br/> CTGTCTGAGATTTTATACCCTTTGACCATCATTAATTTTTTTTTTTTTTTTTTTGAGACA<br/> GTCTCACTCTGTACCCAGTT</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Restriction Sites:</b>           | NotI-NotI                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>ACCN:</b>                        | NM_015962                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Insert Size:</b>                 | 2820 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>OTI Disclaimer:</b>              | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Components:</b>                  | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Reconstitution Method:</b>       | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>RefSeq:</b>                      | <u><a href="#">NM_015962.4</a></u> , <u><a href="#">NP_057046.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>RefSeq Size:</b>                 | 2409 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>RefSeq ORF:</b>                  | 597 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

**Locus ID:** 51077

**UniProt ID:** [Q9Y324](#)

**Cytogenetics:** 14q24.3

**Domains:** PINc, DUF652

**Gene Summary:** Essential protein involved in pre-rRNA processing and 40S ribosomal subunit assembly.  
[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) encodes the longer isoform (1).