

Product datasheet for **SC127586**

C15orf17 (BC015637) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C15orf17 (BC015637) Human Untagged Clone
Tag:	Tag Free
Symbol:	C15orf17
Synonyms:	FLJ00005
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for BC015637, the custom clone sequence may differ by one or more nucleotides

ATGGCGACCGCGGAGCCAGCGGGCGCGGTTGCGGTTGTCTACCCCGGGACCCGGCCAGCGGGGCTC
 GGGACCGCGCGCCGGAGCTGCGGGGCCACCTCCGGGCAGATCGGTAATAGAGCCCTCCGTCTGGGGGA
 GCGCACCCCGCGGCGGTGGAAGCGGGGGCCATACATGGTGACGCGCGCACCTCCATTCAAGCCAAG
 CTGCAGAAGCACCGGACCTGGCCAAGGCCGTTCTGCGGAGAAAAGGCATGCTGGGGGCTCGCCGAACC
 GCCCAGACTCTTCAGGGAAAAGGTCAGTGAAGTTTAACAAGGGCTATACTGCTCTTAGCCAGAGTCCAGA
 TGAACCTGGTGTCCCTCGACTCTGACAGTGATGGGGAGCTGGGATCCAGATACTCCTCCGGGTATTCA
 TCTGCAGAGCAGGTGAACCAGGATGTGAGCCGGCAGCTGCTTCAGGATGGGTATCACCTGGATGAGATTC
 CAGATGACGAGGACTTGGACCTCATCCCCCTAAGCCCATGGCCTCTTCAACATGCTCCTGCTGCTGGTG
 CTGCTTGGGGACTCTTCTTCTGTACCTCCAGTAG


[View online »](#)

5' Read Nucleotide Sequence:	>OriGene 5' read for BC015637 unedited TAATACGACTCACTATAGGGCGGCCGCGAATTCGGCACGAGGATCCGGCAGGGCGGGGCC CGGCCATGGCGACCGCGGAGCCAGCGGGCGCGGTTGCGGTTGTCTACCCCGGGACCCC GGCCCAGCGGGGCTCGGGACCGCGCGCGGGAGCTGCGGGGCCACCCTCCGGGCAGATCG GTAATAGAGCCCTCCGCTGCGGGGAGCGCACCCCGCGGCCGTGGAAAAGCGGGGGCCAT ACATGGTGACGCGCACCCCTCCATTCAAGCCAAGCTGCAGAAGCACCGGGACCTGGCCA AGGCCGTTCTGCGGAGAAAAGGCATGCTGGGGGCCTCGCCGAACCGCCAGACTCTTCAG GGAAAAGGTCAGTGAAGTTTAACAAGGGCTATACTGCTCTTAGCCAGAGTCCAGATGAAA ACCTGGTGTCCCTCGACTCTGACAGTGATGGGGAGCTGGGATCCAGATACTCCTCCGGGT ATTCATCTGCAGAGCAGGTGAACCAGGATGTGAGCCGGCAGCTGCTTCAGGATGGGTATC ACCTGGATGAGATTCCAGATGACGAGGACTTGGACCTCATCCCCGCTAAGCCCATGGCCT CTTCAACATGCTCCTGCTGCTGGTGCTGTCTTGGGACTCTTCTCCTGTACCCTCCAGT AGACATTACCCTCAAGGTGGGCCCTGCTCACCATGATGCCTACAGAGCCCTGTTGGCCA CTCATAGCTCACACAGTGCAGTGAGCACTGAAGTCACTGACCCCTAGAAAGTGACCTGCA AGTAGCCCTTACCNTGTGCTCTCAGCCCAGAGCCCCATCTGGGCTGTGGGCTGAAGGTA NAATNCTGTANCTCTCTGCTGTGNAGCANAGTTGAAAANTCACTTATAGGGCAGACAA NATCTGGTCTATGGTGGTGACAAGTGCTAGGAACCTT
Restriction Sites:	NotI-NotI
ACCN:	BC015637
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	BC015637.2 , AAH15637.1
RefSeq Size:	2322 bp
RefSeq ORF:	597 bp
Locus ID:	57184
Cytogenetics:	15q24.1-q24.2
Protein Families:	Druggable Genome