

Product datasheet for **SC200005**

GADD45GIP1 (NM_052850) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	GADD45GIP1
Synonyms:	CKBBP2; CKbetaBP2; CRIF1; MRP-L59; PLINP; PLINP-1; Plinp1; PRG6
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_052850
Insert Size:	1119 bp
Insert Sequence:	>SC200005 3'UTR clone of NM_052850 The sequence shown below is from the reference sequence of NM_052850. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CCAGCAGCCTCTGGGGCACCCAGCTCCAGGCTTTGTCCCTTCCCAATAAAGCCTGCTACCTGGCAGT
ACCCCTGAAGAGATCTGCGTCTCTCAAGGTGTGCCTCTCTGGGTCCTTGGTTCCTCCCCAATACTT
AAGAGCTGGGGAATCCCCCACTTGGCTGTACTCAAAAGGGGACAGGGACTCTAGGCCCTGCTGACTAAC
TGGGGCTTCAGTGGGAGGAAAAGGGGGGTGTACAATGGGCAAATAATGGAGAAGCAGATGGCAGAAGAG
AAGAAGGCTTCTGTTTACCAACATTCTCAGTAATTAGCCAATTACAAGGCTGGTGGGGGGGGAAG
TATAGCCAAACAGACTGTGGTGATGATGGTGGGGCCGGGAGAGCAGTCCTTGTAACGTTGGGCTG
GGGTCCACTGCCTGGAATGACACACTGGGGGAGGAGACCCGGAAGCCTCAATCTGCAGGGACTCACA
GACTCCTCTCCCCCTTCCAGGCACCCGTCACCATGCCTGCCTGTGCCAGCCAGCCCCCAAGACCCAG
ACAGGTGGGCGCGGGGAGCCCCAGCAGGCACAGACCATAGGGCTAGTTTTTAAACTTTATTCACTTCA
AAACCTTTATCAGAGACACGGTTCTGTTCTGGGGTGGGGGTGGCCTTGACAATGAGCTGATGTTAGCCC
TTCTCCAGGCTGGAGGCCAAGGCCAGACTGGGGAGGGGGCCTTGGGCCCTTCAGGACTTGCCTAGATG
GAAAGATTTGGGCCACTCCAATGGAGGATCAAAGTTGGGGCTCCAGCTTCCCCTCATCTAGAAAGGGA
TGGGGGGGGTCTCAGGGGTTTTTCGGAGGATGGAGAGGCAAGCAGCCGCTTCCGTCTGGCTGCAG
AGAGCTGGGACAGAGGCCAAGAGGGGCCCATCTGTGCCTCCAGTCCTGCCAGCTTCTCTGAGCAGGC
TGGGCTGGGAGGAGGACGGGGGACAGCTGTCTGGCCGGGAGATGGGAACAGGGGCCACTTATAGGAG
GACAGGAAGCTGGACTGTTTGCTTGTTTTTTAAGATTTTGTATTTTTTCCATTTCTAAATAAAC
ATGAATATATATATA
ACGCGTAAGCGGCCCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG

```



Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_052850.4
Summary:	This gene encodes a nuclear-localized protein that may be induced by p53 and regulates the cell cycle by inhibiting G1 to S phase progression. The encoded protein may interact with other cell cycle regulators. [provided by RefSeq, Aug 2012]
Locus ID:	90480
MW:	40.4