

## Product datasheet for **SC200026**

### **RPL6 (NM\_001024662) Human 3' UTR Clone**

#### Product data:

Product Type:	3' UTR Clones
Product Name:	RPL6 (NM_001024662) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	RPL6
Synonyms:	L6; SHUJUN-2; TAXREB107; TXREB1
ACCN:	NM_001024662
Insert Size:	65 bp
Insert Sequence:	>SC200026 3'UTR clone of NM_001024662 The sequence shown below is from the reference sequence of NM_001024662. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGAATTTATCCTCACAAATTGGTGTTCBAATGTCTTAAGAACCTAATTAATAGCTGACTACAT ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 µg dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001024662.3</a></u>



[View online »](#)

**Summary:** This gene encodes a protein component of the 60S ribosomal subunit. This protein can bind specifically to domain C of the tax-responsive enhancer element of human T-cell leukemia virus type 1, and may participate in tax-mediated transactivation of transcription. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed throughout the genome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]

**Locus ID:** 6128

**MW:** 2.5