

Product datasheet for **SC200063**

SLC41A3 (NM_001008486) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SLC41A3 (NM_001008486) Human 3' UTR Clone
Symbol:	SLC41A3
Synonyms:	SLC41A1-L2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001008486
Insert Size:	80 bp
Insert Sequence:	<p>>SC200063 3'UTR clone of NM_001008486</p> <p>The sequence shown below is from the reference sequence of NM_001008486. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <p>GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCCGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TTTTGTGTGTAGCAAAACCTCTATTAGATTTCATAATCAGAGAAGTGTAATAAACAGATTAT ATTGTACTTGA ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG</p>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001008486.2</u>
Locus ID:	54946
MW:	3


[View online »](#)