

Product datasheet for **SC200242**

ODF3B (NM_001014440) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ODF3B (NM_001014440) Human 3' UTR Clone
Symbol:	ODF3B
Synonyms:	FAP123; ODF3L3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001014440
Insert Size:	84 bp
Insert Sequence:	<p>>SC200242 3'UTR clone of NM_001014440</p> <p>The sequence shown below is from the reference sequence of NM_001014440. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <p>GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC GCCCGCTGGTGACCGACGCGGACAACTGACCCGCCAGGCGGGAGCGGCCCCACAGTGTTTGCTTAAA GTCTGCGAGTCCGCA ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG</p>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001014440.4</u>
Locus ID:	440836
MW:	3.1


[View online »](#)