

## Product datasheet for **SC200306**

### **RPS21 (NM\_001024) Human 3' UTR Clone**

#### Product data:

Product Type:	3' UTR Clones
Product Name:	RPS21 (NM_001024) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	RPS21
Synonyms:	HLDF; S21
ACCN:	NM_001024
Insert Size:	84 bp
Insert Sequence:	>SC200306 3'UTR clone of NM_001024 The sequence shown below is from the reference sequence of NM_001024. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GCCGATGGCATCGTCTCAAAGAAGCTTTGACTGGAGAGAATCACAGATGTGGAATATTTGCATAAATA AATAATGAAAACCTA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001024.4</a></u>



[View online »](#)

**Summary:** Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S21E family of ribosomal proteins. It is located in the cytoplasm. Alternative splice variants that encode different protein isoforms have been described, but their existence has not been verified. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]

**Locus ID:** 6227

**MW:** 3.2