

Product datasheet for **SC200339**

MYH (MUTYH) (NM_001048174) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MYH (MUTYH) (NM_001048174) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	MUTYH
Synonyms:	MYH
ACCN:	NM_001048174
Insert Size:	94 bp
Insert Sequence:	>SC200339 3'UTR clone of NM_001048174 The sequence shown below is from the reference sequence of NM_001048174. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GCACACAGCCTCAACAGTGCAGCCAGTGCACCTCTGAAAGCCCCATTCCCTGAGAATCCTGTTGTT AGTAAAGTGCTTATTTTTGTAGTTA ACGCGT AAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001048174.2</u>



[View online »](#)

Summary: This gene encodes a DNA glycosylase involved in oxidative DNA damage repair. The enzyme excises adenine bases from the DNA backbone at sites where adenine is inappropriately paired with guanine, cytosine, or 8-oxo-7,8-dihydroguanine, a major oxidatively damaged DNA lesion. The protein is localized to the nucleus and mitochondria. This gene product is thought to play a role in signaling apoptosis by the introduction of single-strand breaks following oxidative damage. Mutations in this gene result in heritable predisposition to colorectal cancer, termed MUTYH-associated polyposis (MAP). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2017]

Locus ID: 4595

MW: 3.4