

Product datasheet for **SC200350**

RPS16 (NM_001020) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RPS16 (NM_001020) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	RPS16
Synonyms:	S16
ACCN:	NM_001020
Insert Size:	166 bp
Insert Sequence:	>SC200350 3'UTR clone of NM_001020 The sequence shown below is from the reference sequence of NM_001020. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CGCGCTCGCTACCAGAAATCCTACCGATAGCCCATCGTGACTCAAACTCACTTGTATAATAAACAGT TTTTGAGGGATTTAAAGTTTCAAGAAGTGTGTGGCCTTATGTGTTGGTATTGGATGTTAAACCAGA CAGAAATCAGTAAACATCCTGGACCTAT ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001020.6</u>



[View online »](#)

Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S9P family of ribosomal proteins. It is located in the cytoplasm. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]

Locus ID: 6217

MW: 6.3