

Product datasheet for **SC200388**

WDR62 (NM_001083961) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	WDR62 (NM_001083961) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	WDR62
Synonyms:	C19orf14; MCPH2
ACCN:	NM_001083961
Insert Size:	110 bp
Insert Sequence:	>SC200388 3'UTR clone of NM_001083961 The sequence shown below is from the reference sequence of NM_001083961. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GCCGTGCGGAGGAAGGCACGGGGGCAC TGA GGGCGCAGCCCTCCACCGCAGCCCTGCTGCTTCTGAGG ACTTAGGTATTTTAAGCGAATAAACTGACAGCTTTGAGGAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001083961.2</u>
Summary:	This gene is proposed to play a role in cerebral cortical development. Mutations in this gene have been associated with microencephaly, cortical malformations, and cognitive disability. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2011]
Locus ID:	284403



[View online »](#)

MW:

4