

Product datasheet for **SC200411**

MYH13 (NM_003802) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MYH13 (NM_003802) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	MYH13
Synonyms:	MyHC-eo; MyHC-IIL
ACCN:	NM_003802
Insert Size:	116 bp
Insert Sequence:	>SC200411 3'UTR clone of NM_003802 The sequence shown below is from the reference sequence of NM_003802. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GACGTGGGCAGCCAGAAGATGGAAGAAATGAGGCTCACCTGATGCTCGTTGCCATGGGACACCTCCGAGA GAGTGGAGGGAAAATGTGTGAGAAATAAATTCTCCTAAATACTCGGA ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_003802.3
Summary:	Fast twitching myosin mediating the high-velocity and low-tension contractions of specific striated muscles.[UniProtKB/Swiss-Prot Function]
Locus ID:	8735
MW:	4.3



[View online »](#)