

Product datasheet for **SC200423**

ZDHH4 (NM_001134388) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZDHH4 (NM_001134388) Human 3' UTR Clone
Symbol:	ZDHH4
Synonyms:	ZNF374
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001134388
Insert Size:	494 bp
Insert Sequence:	>SC200423 3'UTR clone of NM_001134388 The sequence shown below is from the reference sequence of NM_001134388. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CCATGTCATGAGAGGAAGAAACAAGAATGACAAGTGTATGACTGCCTTTGAGCTGTAGTTCCCGTTTAT
TTACACATGTGGATCCTCGTTTTCCAAGCATGGCTTGTGTTTGTGTTTCTGCTGTGCTTATAAATCA
CTTTCGGTGGCAAGGGAGAGAGGGGAAAATGGGTGTTGACTGAGGAATCCCCCTTGCTTGTCTTCTTT
TGAAACCGGGCATCTCTGAAGTCTGGTGTCAAGGGGATCAAGAGATGACTTCTCAGAGGTTCTAGGTG
ATGCTGAGACCTTGGTGTCTCTAAACTCTGGGCATGTGGACAGGAGGGGCTTGCGGCCGTGCTCTGAC
CTGTGTGATGTGACAGGAGGTCTCATTGACTCAGCGCCTGCGCGTTGTGCTGGCTGTGCTCTCTTTTA
TGCCCCCTCTATTCTCTCTCTCCCCAGGGGATTTTCATCTCAACAACAGAGTGTCTCACCCAGGGC
TGCCTTCGGC
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001134388.2
Summary:	Palmitoyltransferase that could catalyze the addition of palmitate onto protein substrates including the D(2) dopamine receptor DRD2.[UniProtKB/Swiss-Prot Function]
Locus ID:	55146
MW:	18.6