

Product datasheet for **SC200433**

TMEM171 (NM_001161342) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TMEM171 (NM_001161342) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	TMEM171
Synonyms:	PRP2
ACCN:	NM_001161342
Insert Size:	109 bp
Insert Sequence:	>SC200433 3'UTR clone of NM_001161342 The sequence shown below is from the reference sequence of NM_001161342. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CCTCTATCTTCTGAGCCTCCCCACCGTAACTATGGACTCTAGTTCAGTTTTATATGCAATGGATCAC TATTTTATTTAATTTTTTTTAAATAAAAAATACAATAGCA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001161342.3</u>
Locus ID:	134285
MW:	4.2



[View online »](#)