

Product datasheet for **SC200485**

IRF7 (NM_004029) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	IRF7 (NM_004029) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	IRF7
Synonyms:	IMD39; IRF-7; IRF-7H; IRF7A; IRF7B; IRF7C; IRF7H
ACCN:	NM_004029
Insert Size:	120 bp
Insert Sequence:	>SC200485 3'UTR clone of NM_004029 The sequence shown below is from the reference sequence of NM_004029. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TTCCTTATGGAGCTGGAGCAGCCCGCCTAGAACCCAGTCTAATGAGAACTCCAGAAAGCTGGAGCAGCC CACCTAGAGCTGGCCCGCGCCCGCCAGTCTAATAAAAAGAACTCCAGAACA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_004029.4</u>



[View online »](#)

Summary: IRF7 encodes interferon regulatory factor 7, a member of the interferon regulatory transcription factor (IRF) family. IRF7 has been shown to play a role in the transcriptional activation of virus-inducible cellular genes, including interferon beta chain genes. Inducible expression of IRF7 is largely restricted to lymphoid tissue. Multiple IRF7 transcript variants have been identified, although the functional consequences of these have not yet been established. [provided by RefSeq, Jul 2008]

Locus ID: 3665

MW: 4.3