

## Product datasheet for **SC200527**

### TRPT1 (NM\_001160389) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	TRPT1 (NM_001160389) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	TRPT1
ACCN:	NM_001160389
Insert Size:	90 bp
Insert Sequence:	>SC200527 3'UTR clone of NM_001160389 The sequence shown below is from the reference sequence of NM_001160389. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TCCAGAGAAAGGAGGAGGATCCAACAATAAATATTAATTTATAAAAAAGAAATTTAAAAAGTAACAA GAAAGAACTCGTTTGAACCA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<a href="#">NM_001160389.2</a>
Summary:	Catalyzes the last step of tRNA splicing, the transfer of the splice junction 2'-phosphate from ligated tRNA to NAD to produce ADP-ribose 1''-2'' cyclic phosphate.[UniProtKB/Swiss-Prot Function]
Locus ID:	83707
MW:	3.9



[View online >](#)