

Product datasheet for **SC200575**

LENG1 (NM_024316) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: LENG1 (NM_024316) Human 3' UTR Clone
Symbol: LENG1
Mammalian Cell Selection: Neomycin
Vector: pMirTarget (PS100062)
ACCN: NM_024316
Insert Size: 590 bp
Insert Sequence: >SC200575 3'UTR clone of NM_024316
The sequence shown below is from the reference sequence of NM_024316. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CGCCAGCAGGACCTCACCTTACTACTTGACTCTCTGAGGGGTACAGGAGAGGCCGCTGCTGCCAGCCG
TCATATAAACTATTTATTCATAAATATTTCCAAAATGAAAATAGGTTTACAAAAATGTCCTCAC
TGGGAGGGGAGGAGGGGGCAGCCCTCGCCCCGGGCCCCAGGGTGGGGCTGAGAGGAAAACCTCCCG
GCCCTCCCTGCTTCTGGGAGAGGGGATGCCCGTGGCTTGGGCCCTCCCTCCAGTCTTCCAGGGC
AGGGCCCTCACCTGGGCAGGGGATCAGCATGCGGGGAAGGGGTGGGTAGAGGGAGGGCCGGTGTCT
ACTGGAGGTCCCGTCTCCAGGTAGCGGTACTCAAAGGTGAAGCCTTCTTCTCCGCTGGCCCCACT
TCTCGTAGTCAAAGTAGATGTAGGTGCCCTGGCCGGGGAGAAGCGGTAGTGTGACGAGGAGGT
GGTCTGGGATCTGGCCGGACCAACAGACAAAGGGGACAATTCTTAGGGCTGTGGATGTGTAGGCACC
GGGCCAGCTGCCCTGCACGCACACTCTCATCCATCC
ACGCGTAAGCGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI
OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_024316.3](#)

Locus ID: 79165

MW: 20.2