

Product datasheet for **SC200585**

MRPL55 (NM_181465) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MRPL55 (NM_181465) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	MRPL55
Synonyms:	AAVG5835; L55nt; MRP-L55; PRO19675
ACCN:	NM_181465
Insert Size:	105 bp
Insert Sequence:	>SC200585 3'UTR clone of NM_181465 The sequence shown below is from the reference sequence of NM_181465. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CGACAGTTCTGGACCAGGACCAAGAAGTGAACCGTGGCTCCAGCCACCCGGGACATTGCTAAGATGGGA GGGCTGTTCTTAAATCACTCGTTCTTGAAGCTGCCA ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_181465.3</u>



[View online »](#)

Summary:

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Multiple transcript variants encoding two different isoforms were identified through sequence analysis. [provided by RefSeq, Jul 2008]

Locus ID:

128308

MW:

4.2