

## Product datasheet for **SC200600**

### **ANKFN1 (NM\_153228) Human 3' UTR Clone**

#### Product data:

Product Type: 3' UTR Clones

Product Name: ANKFN1 (NM\_153228) Human 3' UTR Clone

Vector: pMirTarget (PS100062)

Symbol: ANKFN1

ACCN: NM\_153228

Insert Size: 131 bp

Insert Sequence: >SC200600 3'UTR clone of NM\_153228  
The sequence shown below is from the reference sequence of NM\_153228. The complete sequence of this clone may contain minor differences, such as SNPs.  
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAACTTGGTATAGTAGCTTGTTTCACCTAGAAATATTAACCCAGCCTCCTTATAATAAAATCACAAAGT
TATATCTGTTCCCTTTGTCCCAGTGGAGGGTCAATAAATCACATGATGGCTTTGGCAACAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-Mlul

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM\\_153228.3](#)

Locus ID: 162282

MW: 4.8



[View online »](#)