

Product datasheet for SC200606

GON4L (NM_032292) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GON4L (NM_032292) Human 3' UTR Clone
Symbol:	GON4L
Synonyms:	GON-4; GON4; YARP
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032292
Insert Size:	486 bp
Insert Sequence:	>SC200606 3'UTR clone of NM_032292 The sequence shown below is from the reference sequence of NM_032292. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATGGATAGTCCTAAAATCATTGTATTGATTTGTGAATGCATTATGGGACATGATTGTGGAGTTGAGG
TGAAATGAGATGGAAAGGATGAAATTTTACTTATTATATTAACCTCGTTTACACATTAAATTGAATCAC
CTGACCATACTTTAAATGCCTTAAAAAATAGTAATAGGCTGGGCTTGGTGGCTCACTCCTGTAGTCCCA
GCACTTTGGGAGGCTGAGGTGGGTGGATCACAAGGCCAGGAGTTTCAGACTGGCCTGGCCGACGTGGTA
AAACCCCATCTATATTAATAAATAACAAAAATTAGCCAGGCATGGTGGCGGACACCTGTAATCCCTGCTG
CTTGGGAGGCAGAGGCAGGAGAATTGCTTGAGCCCGAGAGACAGAGGTTGCAGTGAGCCGAGATCACGT
CATTGCTCTCCAGCCTGGGCAACAGAGCAAGGCTCTGTCTCAAAAAAAAAAAAAAAAAAAAAAAAAAAAA
AAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).


[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_032292.6
Summary:	Has transcriptional repressor activity, probably as part of a complex with YY1, SIN3A AND HDAC1. Required for B cell lymphopoiesis.[UniProtKB/Swiss-Prot Function]
Locus ID:	54856
MW:	18.7