

Product datasheet for **SC200616**

FLAD1 (NM_025207) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FLAD1 (NM_025207) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	FLAD1
Synonyms:	FAD1; FADS; LSMFLAD; PP591
ACCN:	NM_025207
Insert Size:	104 bp
Insert Sequence:	>SC200616 3'UTR clone of NM_025207 The sequence shown below is from the reference sequence of NM_025207. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GAAGAAGAGGAGCGGAACCTCCGCACATGACCTCCACCCTAGGAGGGAGGGAAGGACACCGTCCTAGG GTATAACCTGGCAATAAACCGTGCCTCTCACTGTG ACGCGT AAGCGGCCGCGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_025207.5
Summary:	This gene encodes the enzyme that catalyzes adenylation of flavin mononucleotide (FMN) to form flavin adenine dinucleotide (FAD) coenzyme. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]
Locus ID:	80308



[View online »](#)

MW: 3.8