

Product datasheet for **SC200629**

IQCN (NM_001145305) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	IQCN (NM_001145305) Human 3' UTR Clone
Symbol:	IQCN
Synonyms:	KIAA1683
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001145305
Insert Size:	112 bp
Insert Sequence:	>SC200629 3'UTR clone of NM_001145305 The sequence shown below is from the reference sequence of NM_001145305. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC TCACGGCACATGCATTGGCCTGGCATC TAG GACCCTGGCTCCCTGCAGTGGGGACTTCGTGGGAGGCAC TCATGGCTCTCTGGGTCTAATGAATAAAGTCTCCACAGCCTA ACGCGT AAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001145305.2
Locus ID:	80726
MW:	4.2



[View online »](#)