

## Product datasheet for **SC200636**

### **C6orf48 (NM\_001040437) Human 3' UTR Clone**

#### Product data:

Product Type:	3' UTR Clones
Product Name:	C6orf48 (NM_001040437) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	C6orf48
Synonyms:	D6S57; G8
ACCN:	NM_001040437
Insert Size:	133 bp
Insert Sequence:	>SC200636 3'UTR clone of NM_001040437 The sequence shown below is from the reference sequence of NM_001040437. The complete sequence of this clone may contain minor differences, such as SNPs. <b>Blue</b> =Stop Codon <b>Red</b> =Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CTTATTGGGGGAAAAATGCCTGAAAAT <b>TAG</b> GGGGCACTTCAAGTAGATAGCTTCTATTTCTATATTTG TCTTATATACAAGTATTTGCTTTTATCAAATAATTCCAATAAAGCATTTTAAAGTAAAGAAGA <b>ACGCGT</b> AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001040437.2</a></u>
Locus ID:	50854
MW:	5.1



[View online »](#)