

Product datasheet for **SC200649**

SGT1 (ECD) (NM_007265) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SGT1 (ECD) (NM_007265) Human 3' UTR Clone
Symbol:	SGT1
Synonyms:	GCR2; HSGT1; SGT1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_007265
Insert Size:	945 bp
Insert Sequence:	>SC200649 3'UTR clone of NM_007265 The sequence shown below is from the reference sequence of NM_007265. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
AGACCAACAAGTAAGCCAACAAAAATTAACAGCACATTTAGCTTCTCTTTTTCTTTTAAATAAAT
ATTGAATATGATTCTGTTCAATTTCTAGTTTAGATTACTCATTTAGTAAAAATGTTTCTTCATTACAAGT
TTTACAGATTGACTTCTCCATGAGCCCCAAATGTGGCCCTCTGGCACAAGTCTTTGTGATCTCAGCTG
TGGGTTATACACAGGTTTTGAACTCGGCAGGGTAGAATCTAGCCCTAGACAACCTTGAGATGAAATTC
ACACACTTCATGCTAAAGACATACCCTACCACAATATACCCTTGGAGTAATTCTAAATATTTAAATTTTC
TCCCTTACTGCTATTCATCCTCAGATCAAAATTTCAAGTTGACAAATTTGGATTGAAAACCTCATTCAC
TTTTTCAATTTTTTAATGAGGCATAACATGCATATGGTTAAGTGCACAAATCTTAAGTGTACAGTCAG
TGAATTTTGCATGTATGTATGTACCTGAAGCAAGATATAGAATTTACTTGTAAATAACAGCTAT
ATTGATATGTCACATACCAAAAAATTTACCCTTTAAATTATACAATTCAGTGGTTTTAGTATATTTAC
AAAATTGTGCAACCATCCCTACCATCTAATTTCTAGAATATTAGAAACCCATACTCACTAATAGTCTCT
CATTCCCTTCTACTCCCAGCCCCTAGAAACCACTAGTCTACTTTCTGTCTCTATGGATTTGCCTATTTT
GGATATTTACATCTGATCTGACTTAATTTGTTTTAATTAATTAATTTGGTTTGATTTAATATAACT
GATATATCTTTTGCTAAGATGAAATGTTTATTTCTTAAATGATGCTCATATTCCTCCTAAAATAAATC
AAGTATTATCTATATATCATAAAATTCACCATTAAGCAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_007265.3
Summary:	Regulator of p53/TP53 stability and function. Inhibits MDM2-mediated degradation of p53/TP53 possibly by cooperating in part with TXNIP (PubMed:16849563, PubMed:23880345). May be involved transcriptional regulation. In vitro has intrinsic transactivation activity enhanced by EP300. May be a transcriptional activator required for the expression of glycolytic genes (PubMed:19919181, PubMed:9928932). Involved in regulation of cell cycle progression. Proposed to disrupt Rb-E2F binding leading to transcriptional activation of E2F proteins (PubMed:19640839). The cell cycle -regulating function may depend on its RUVBL1-mediated association with the R2TP complex (PubMed:26711270). May play a role in regulation of pre-mRNA splicing (PubMed:24722212).[UniProtKB/Swiss-Prot Function]
Locus ID:	11319
MW:	37.4