

Product datasheet for **SC200651**

C3orf78 (SMIM4) (NM_001124767) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C3orf78 (SMIM4) (NM_001124767) Human 3' UTR Clone
Symbol:	C3orf78
Synonyms:	C3orf78
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001124767
Insert Size:	1435 bp



[View online »](#)

Insert Sequence:	>SC200651 3'UTR clone of NM_001124767 The sequence shown below is from the reference sequence of NM_001124767. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CAGTATCAGAGAAGGCTGGAAGATGAATGAGACTGAACTTCAGCAGTCAATAAAGTCAATATGAATTTT TACTATTGGTTTTAGTGCCTTTTAAATATACTTTTCAGATCTCTTGTCTGACACAGCTTGTCTTTA TCAGCTGCACAATATTCCAATACTGTTTTGCAAAGAATTTTATAGTTCTGTCTCTGCTTGAGTACT TCCCAAGTCACTGCTGGACAGCACTTGTGGGAAAATTATTTCTCCTAAAACACTCACTCGATGAATCT AGTCCATGAAATATTTCTTAAACAATTCAAGGATCAATCAGGATGCTTTTCAGTGAAGCCGAATTGTT TTCGGCAATCTTTGTTTTAGGCAACTAATTATAACCTTATTTACTTTTTTGAGAATTCAAATTAAT TATGGACACATTCTCTATTGCCAGATCATTTAAACCAGGATAATTATGCAGGTTGATTTCTGCCAGCT ACTAGACTATAAACTTTTTGCTAGGGTGGTGGAAATCTTGATTGACTACGCATTTCTCATGCTAGTGT TCTTGTGAACAGAAATATGCTCCTCCCTTGCTTGCTCCCTAAACTGCAAGGGCTGGATCTTCATGAG GACTGTGGGCAACAGGTTTCACACTACTCTAGGATGAAGATGTTGCTGTTTCGATTGAGTCTAGAAC AGTTGCTCCTTCCATGAGAAAGTGAACCTACACAGTATAAAGCATCCCAGACTTGAGCTAAATCAGAGG CTCAGTCTATTGAGTGATTTGGCAATCTGACCCAGTGTAAGTTTGGATCAGCCTGGTAGAAAAAGCAGA TACCTTTCCATGCTTGCATAGCCACATCTAACATCGTTTCAAGCAAAAAAAAAATCTACATACTAAAGAC TAAGTTTCAGAATAGTTTCTAACTATTGCATGAGAATAAAGCTGATGCTTAAGGCCATAGGGTATCCTC AAGTTAATGATCTCTGCCACAAACATGTCTCCTAATTTTTCTAGATTTTTCTTCTAATTTTAGTCCAC TTATCTTAGAACCAATTGAGTAAAGAACCTAACTAACTACTCCATCTCCTTAAAGCACTACCTACA CTACAGTGCAGCAAGCCAGAGTCACAGCTGGGAGATGGCTTGAATGTGAAAGGCCAGTAATGACCTCT GCCATCAGGGAGAGGACGGTGGTTCATCGTATTCTCATGCCAGTTGTGTCATCTAGGATTCTGGTAA ACCTAACCTGATAAACCATTTTTAGTGACTTGACTGCTGTATAAAAACTATAGTATGATTCTCAATTAT ATGGCCATACTCCATGTTTGCCTAACCACCTTGTGTTAGGAGTTAGTATTTGATACATATAGAAAA TACGTCATACCTTTTACAAAGAATAATAGTAAACACTAATAAAACCATTGGCTAA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001124767.2</u>
Locus ID:	440957
MW:	55.1