

Product datasheet for **SC200706**

ADRM1 (NM_175573) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ADRM1 (NM_175573) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ADRM1
Synonyms:	ARM-1; ARM1; GP110; PSMD16
ACCN:	NM_175573
Insert Size:	147 bp
Insert Sequence:	>SC200706 3'UTR clone of NM_175573 The sequence shown below is from the reference sequence of NM_175573. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GACGAAGAGGAGGACATGAGCCTGGACTGAGCCACGCGCCGTCTCCGAGGAACTGGGCGCTTGCAGTG CGTTGCACACCCTCACCTCCCACCCACTGATTATTAATAAAGTCTTTTCTTTTACCTGCCAAAAAAAA AAAAAAAA ACGCGT AAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_175573.2</u>



[View online »](#)

Summary: This gene encodes a member of the adhesion regulating molecule 1 protein family. The encoded protein is a component of the proteasome where it acts as a ubiquitin receptor and recruits the deubiquitinating enzyme, ubiquitin carboxyl-terminal hydrolase L5. Increased levels of the encoded protein are associated with increased cell adhesion, which is likely an indirect effect of this intracellular protein. Dysregulation of this gene has been implicated in carcinogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]

Locus ID: 11047

MW: 5.7