

## Product datasheet for **SC200729**

### STK32A (NM\_145001) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	STK32A (NM_145001) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	STK32A
Synonyms:	YANK1
ACCN:	NM_145001
Insert Size:	137 bp
Insert Sequence:	>SC200729 3'UTR clone of NM_145001 The sequence shown below is from the reference sequence of NM_145001. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GATACCTGGCTCTCCTACAAGTCCCACATGATTGGAGTTTCAGGAGACCGAAGCCCAGGCACATGTATT TTGCAAACTACTGAAGTTTCTGATAATGACGGATATCAACAATTAACGCTTACTTCTTGCCAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<a href="#">NM_145001.4</a>
Locus ID:	202374
MW:	5.1



[View online »](#)