

Product datasheet for **SC200742**

UBAPIL (NM_001163692) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones

Symbol: UBAPIL

Synonyms: UBAP-IL

Mammalian Cell Neomycin

Selection:

Vector: pMirTarget (PSI00062)

ACCN: NM_001163692

Insert Size: 367 bp

Insert Sequence: >SC200742 3'UTR clone of NM_001163692
The sequence shown below is from the reference sequence of NM_001163692. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
CTGGAGGAGCTGGTGGCCTGTGCCAGAGACCAGGAGGCACTCCACGATGCCTGGGCATCTAGGCCCA
ATGCCAGACCTGGCAACCCATCCTAACTGTATTATTAACCAAGTAAGAGTAACAAAGTGATGCATAC
TTATTGTAAAAACAAATAATACCAAGGGTGACAGTTAACAGCAGCAGAATCTGCTGAAAGAGTTGAA
ATAGTTGCCTCTGCAGAGGGGAGATGGGAGAGCCGGACGGGGAGCTATGTTCTTCATGCATGCCTGTG
AAACCACCTGATTGACTCTTTACACTGTGTGCATTAGTAACCTTGATTTTAAGAAAATTAACAAAT
AAATAATGGGCAGTTCTTCTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001163692.2</u>
Locus ID:	390595
MW:	13.5