

Product datasheet for **SC200776**

C16orf13 (METTL26) (NM_001040160) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C16orf13 (METTL26) (NM_001040160) Human 3' UTR Clone
Symbol:	C16orf13
Synonyms:	C16orf13; JFP2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001040160
Insert Size:	126 bp
Insert Sequence:	>SC200776 3'UTR clone of NM_001040160 The sequence shown below is from the reference sequence of NM_001040160. The complete sequence of this clone may contain minor differences, such as SNPs. Blue =Stop Codon Red =Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA GCGATCGCC CTCTGTGAGGCACGAACCCTGCCTCCCT AG GCCGGACCTTGTGGACGACAGCCCCACCCAGTCTGTGCT CTCAGCCGCTGGCCGAAGGGCCAGCCTGCTCAGAATAAAGCATGTCCTGCTGCCGG ACGCGT AAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001040160.3
Locus ID:	84326
MW:	4.5



[View online »](#)