

Product datasheet for **SC200779**

DNA helicase HEL308 (HELQ) (NM_133636) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	DNA helicase HEL308
Synonyms:	HEL308
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_133636
Insert Size:	137 bp
Insert Sequence:	>SC200779 3'UTR clone of NM_133636 The sequence shown below is from the reference sequence of NM_133636. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GCTGTGGCTTCTTCCACTGACAAAGCATGAAGCTATCTGATGAGAATTTTCATATATGAATTATTTATT GTACATGTGTGTTATTAAGATAACTTATACTTTATAAATGAAAAATACATTTTCAGTTATACATTAA ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.



RefSeq: NM_133636.5

Summary: HEL308 is a single-stranded DNA-dependent ATPase and DNA helicase (Marini and Wood, 2002 [PubMed 11751861]).[supplied by OMIM, Mar 2008]

Locus ID: 113510

MW: 5.4