

Product datasheet for **SC200801**

DNAJC4 (NM_005528) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	DNAJC4
Synonyms:	DANJC4; HSPF2; MCG18
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_005528
Insert Size:	127 bp
Insert Sequence:	>SC200801 3'UTR clone of NM_005528 The sequence shown below is from the reference sequence of NM_005528. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GAGATCGTGCCCCGGGGCGCCGGCCCCGAGGGGCTCACCTGGATGGGGCTGCAGTGC GTTCCCGCTT TGCTTCCTTCCCTGGACGGCCCGCTCCCCGAAACGCGCGCAATAAAGTGATTTCGCAGA ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.



RefSeq: NM_005528.4

Locus ID: 3338

MW: 4.4