

## Product datasheet for SC200809

### Lipolysis Stimulated Lipoprotein Receptor (LSR) (NM\_205835) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	Lipolysis Stimulated Lipoprotein Receptor (LSR) (NM_205835) Human 3' UTR Clone
Symbol:	Lipolysis Stimulated Lipoprotein Receptor
Synonyms:	ILDR3; LISCH7
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_205835
Insert Size:	135 bp
Insert Sequence:	<p>&gt;SC200809 3'UTR clone of NM_205835</p> <p>The sequence shown below is from the reference sequence of NM_205835. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> <pre> GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA<del>CGATCGCC</del> GCCCTGAGTCGGGAAAGTTTAGTCGTC<del>TGA</del>TCTGACGTTTTCTACGTAGCTTTTGTATTTTTTTTTTTA ATTTGAAGGAACACTGATGAAGCCCTGCCATACCCCTCCCGAGTCTAATAAAACGTATAATCACAA <del>ACGCGT</del>AAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG           </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<a href="#">NM_205835.4</a>


[View online »](#)

<b>Summary:</b>	Probable role in the clearance of triglyceride-rich lipoprotein from blood. Binds chylomicrons, LDL and VLDL in presence of free fatty acids and allows their subsequent uptake in the cells (By similarity).[UniProtKB/Swiss-Prot Function]
<b>Locus ID:</b>	51599
<b>MW:</b>	5.2