

## Product datasheet for **SC200851**

### CCDC169 (NM\_001144981) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CCDC169 (NM_001144981) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CCDC169
Synonyms:	C13orf38
ACCN:	NM_001144981
Insert Size:	448 bp
Insert Sequence:	>SC200851 3'UTR clone of NM_001144981 The sequence shown below is from the reference sequence of NM_001144981. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site  GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAA <b>ACGATCGCC</b> CCAAACCATCTCCAGAACTCCATCCA <b>TG</b> ATTCCAGAACTGGATGGATTCTATTTCTTCTTTTTTTT TCCATTTATAATATTCAAGTTATAAAGTTAATGGTCAGCACATTTCTTCAGGAATAAAATAGATAAGA ACTGCCACATTTTTCTTTAGTTCCTTTTCCAGAGGTGTGTTCTTTGTCTGCATAAAAAGTAATTTCTC CAGTGACCCTAAACTATGTAATAGAAGTCTTGAGAGTTTTAAAAATATCCTAGGGAATCTGAAGGCTGG CATCATTACACTTGCTTTCAATAATATACCCATGATCTGCTCTGATAGTAACGTTGCTATTCTATTTGG TTTCTAGTTGACTGTTAGCAAAAAGTGTGATTTCATAAAACCTTTTTGCCATTTTGGGCTGGATCTTGC TCAGTTTTTTAAATTGCAGGTTTAAACAAAGCCG <b>ACGCGT</b> AAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u><a href="#">NM_001144981.3</a></u>



[View online »](#)

Locus ID: 728591

MW: 17.9