

Product datasheet for **SC200852**

CCDC169 (NM_001144986) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CCDC169 (NM_001144986) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: CCDC169
Synonyms: C13orf38
ACCN: NM_001144986
Insert Size: 448 bp
Insert Sequence: >SC200852 3'UTR clone of NM_001144986

The sequence shown below is from the reference sequence of NM_001144986. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CCAAACCATCTTCCAGAACTCCATCCAATGATTCCAGAACTGGATGGATTCTATTTCTTCTTTTTTTT
TCCATTTATAATATTCAAGTTATAAAGTTAATGGTCAGCACATTTCTTCAGGAATAAAATAGATAAGA
ACTGCCACATTTTCTTTAGTTCCTTTTCCAGAGGTGTGTTCTTTGTCTGCATAAAAAGTAATTTTCTC
CAGTGACCCTAAACTATGTAATAGAAGTGCTTGAGAGTTTTAAAAATATCCTAGGGAATCTGAAGGCTGG
CATCATTACACTTGCTTTCAATAATATACCCATGATCTGCTCTGATAGTAACGTTGCTATTCTATTTGG
TTTCTAGTTGACTGTTAGCAAAAAGTGTGATTTCATAAAACCTTTTTGCCATTTTGGGCTGGATCTTGC
TCAGTTTTTTAAATTGCAGGTTTAAACAAAGCCG
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001144986.3](#)



[View online »](#)

Locus ID: 728591

MW: 17.9